

evercare[®] Paediatric absorbent drape, 70 x 60 cm

REF: 1463-50

evercare[®] Paediatric adhesive absorbent drape, sterile is a versatile, impervious, extra absorbent drape. It is especially designed for patients with sensitive skin as the adhesive edge is only 2,5 cm wide and the adhesive more gentle to the skin. The fabric is also thinner, yet highly absorbent, and conforms nicely to the body contours. Patients that can especially contribute positively from this features are elderly or neonatal.

Key Features

- ▶ Extra absorbent
- ▶ Softer adhesive allows for gentler removal and repositioning of the drape
- ▶ Complies to EN 13795-1 High performance less critical product area requirements




70 x 60 cm

REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1463-50	Paediatric adhesive absorbent drape, sterile, Blue	60 x 70 cm, Adhesive edge 2.5 cm	40 / 160

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-01
EN ISO 10993-5
EN ISO 10993-7
EN ISO 10993-10

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR class Is

Compliance with product standards: EN13795-1

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Laminate	PP nonwoven, blue / PE film, blue

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PE/PET film	250 x 390
Dispenser Pack	TCF or ECF quality	285 x 242 x 287
Transport Carton	TCF or ECF quality	500 x 310 x 590
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1463-50	06438129118704	06438129318708	06438129518702

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.