


# evercare<sup>®</sup> Suture removal set, nonwoven, sterile

REF: 1892-01

evercare<sup>®</sup> Suture removal set, nonwoven, sterile is a compact disposable sterile set offering convenience containing items needed for suture removal.



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
1892-01	Suture removal set, nonwoven, sterile	N/A	15 / 240

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-05 EN ISO 10993-07 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class Is
<b>Compliance with product standards:</b>	EN ISO 11135
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
1 Forceps, 127 mm	Polystyrene, green
3 Triangular swabs, 35 x 35 x 35 mm	Nonwoven, PES/viscose
1 Suture removal knife, 110 mm packed in peel-open plastic/ aluminium bag	Carbon steel
1 Tray	Polystyrene, semitransparent

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	3 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, Plastic film	195 x 80
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	280 x 180 x 100
<b>Transport Carton</b>	Carton, TCF or ECF quality	590 x 390 x 430
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>1892-01</b>	6438129121469	6438129321463	6438129521467

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.