

3-Part Syringes / Regular



Our Terumo three-part syringes are of the highest quality, from development to packaging.

Their highly transparent barrel guarantees easy recognition of air bubbles and the clear printing allow for precise dosage management. The syringes contain no natural rubber and are latex-free. They are environmental-friendly and are designed to have a minimum of waste.

Our range includes Luer slip tip (concentric and eccentric), Luer Lock tip and Catheter tip.

3-Part Syringes / Regular Product Range

Without needle

Reference	Description ml	Type	Unit Box	Shipping Carton	Shelf Life (years)	Sterilization
BS-01T	1	Slip Tip	100	1800	5	EtO
SS-02S	2.5	Slip Tip	100	1200	5	EtO
SS*02LE1	2	Luer Lock	100	2400	5	e-beam
SS+03L1	3	Luer Lock	100	1800	5	e-beam
SS*05SE1	5	Slip Tip	100	1000	5	EtO
SS*05LE1	5	Luer Lock	100	1600	5	e-beam
SS-10ES	10	Eccentric	100	500	5	EtO
SS*10LE1	10	Luer Lock	100	800	5	e-beam
SS+20ES1	20	Eccentric	50	600	5	e-Beam
SS+20L1	20	Luer Lock	50	600	5	e-Beam
SS*30ESE1	30	Eccentric	50	400	5	e-beam
SS*30LE1	30	Luer Lock	50	400	5	e-beam
BS-50C	50	Catheter Tip	25	100	5	EtO
BS-50ES	50	Eccentric	25	100	5	EtO
BS-50LB	50	Luer Lock Barrel Opaque	25	100	5	EtO
BS-50LG	50	Luer Lock	25	100	5	EtO