

Overall with hood
REF 210810, 210811, 210812, 210813

Version: 18.12.2013, rev 1

Printed 2.9.2016

REF	description	size																																	
210810	overall with hood	M																																	
210811	overall with hood	L																																	
210812	overall with hood	XL																																	
210813	overall with hood	XXL																																	
<p>[NATURAL LATEX FREE]</p>																																			
<p>Raw material Main material Sewing thread Zipper Elastic at ankles, cuffs & waist</p>			<p>Spunbond PP nonwoven Nylon</p>																																
<p>REACH regulations:</p>			<p>The products covered by this data sheet do not contain carcinogenic, mutagenic or toxic substances, including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).</p>																																
<p>Features:</p>			<p>A lightweight nonwoven overall for various uses in healthcare settings</p>																																
<p>Shelf life:</p>			<p>3 years</p>																																
<p>Pack counts & carton sizes</p>			<table border="1"> <thead> <tr> <th>REF</th> <th>size</th> <th>pcs/unit</th> <th>pcs/transp.</th> <th>meas. transp.</th> </tr> <tr> <td></td> <td></td> <th>pack</th> <th>ctn</th> <th>ctn mm</th> </tr> </thead> <tbody> <tr> <td>210810</td> <td>M</td> <td>1</td> <td>50</td> <td>560x320x370</td> </tr> <tr> <td>210811</td> <td>L</td> <td>1</td> <td>50</td> <td>560x320x370</td> </tr> <tr> <td>210812</td> <td>XL</td> <td>1</td> <td>50</td> <td>560x320x370</td> </tr> <tr> <td>210813</td> <td>XXL</td> <td>1</td> <td>50</td> <td>560x320x370</td> </tr> </tbody> </table>			REF	size	pcs/unit	pcs/transp.	meas. transp.			pack	ctn	ctn mm	210810	M	1	50	560x320x370	210811	L	1	50	560x320x370	210812	XL	1	50	560x320x370	210813	XXL	1	50	560x320x370
REF	size	pcs/unit	pcs/transp.	meas. transp.																															
		pack	ctn	ctn mm																															
210810	M	1	50	560x320x370																															
210811	L	1	50	560x320x370																															
210812	XL	1	50	560x320x370																															
210813	XXL	1	50	560x320x370																															
<p>Sales unit:</p>			<p>Transport carton</p>																																
<p>Unit pack</p>			<p>PE-bag, 360 x 290 mm, transparent</p>																																

Overall with hood
REF 210810, 210811, 210812, 210813

Version: 18.12.2013, rev 1

Printed 2.9.2016

Transport carton:	Double wall corrugated fibreboard	
Packing tape:	PVC-free	
Compliance with legal requirements:	Personal Protective equipment 89/686/EEC, Cat. I	
Quality standards:	The product is manufactured in compliance with ISO 13485	
Label information:	In accordance with EN ISO 15223-1, EN 980 and EN 1041	
Storage recommendation:	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone	
Disposal recommendation:	Follow local instructions	
Risk if reused:	The product is single use. If reused performance of product may deteriorate, cross contamination may occur	
Barcodes:	Dispenser	Transport carton
N/A	N/A	6438129550368 210810 carton  6 438129 550368
N/A	N/A	6438129553024 210811 carton  6 438129 553024
N/A	N/A	6438129550399 210812 carton  6 438129 550399
N/A	N/A	6438129553451 210813 carton  6 438129 553451