

**Beard cover**  
**REF 210855**

Version: 8.4.2016, rev 1

Printed 8.4.2016

REF	Size	
210855	48 x 24 cm, loop length unstretched 24 cm	
<b>Color:</b> white  <b>[ NATURAL LATEX FREE ]</b>		
<b>Raw material:</b> Body Sewing thread Elastic loop		polypropylene nonwoven, white terylene latex free
<b>REACH regulations:</b>		The product covered by this data sheet does not contain carcinogenic, mutagenic or toxic substances, including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Features:</b>		Non-woven, disposable beard cover. Practical protection for use in various industries for prevention of contamination.
<b>Pack count:</b>		100 pcs/bag, 10 bags/transport carton
<b>Sales unit:</b>		Transport carton
<b>Dispenser:</b>		PE bag 340 x 220 mm
<b>Transport carton:</b>		Cardboard 320 x 240 x 250 mm
<b>Quality standards:</b>		The product is manufactured in compliance with ISO13485
<b>Label information:</b>		In accordance with EN ISO 15223-1, EN 980 and EN 1041
<b>Storage recommendation:</b>		Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone
<b>Disposal recommendation</b>		Follow local regulations

# Product Data Sheet

Protective apparel




Beard cover  
REF 210855

Page 2 (2)

Version: 8.4.2016, rev 1

Printed 8.4.2016

Barcode	Dispenser	Transport carton
210855	N/A	6438129552645 210855 carton  6 438129 552645