



# SELEFA® Shoe covers

REF: 211540, 211541 and 211545

SELEFA® Shoe cover is a practical and economical product for infection control. It is an easy to wear non sterile product that comes in different sizes and colors and thickness to fit all needs. Also suited for non-medical uses.

## Key Features:

- ▶ Easy to put on
- ▶ Different colors and sizes available
- ▶ Single use



REF 211540



REF 211541  
REF 211545

REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
211540	SELEFA® Shoe cover, white, embossed	40 x 15 cm	50/1000
211541	SELEFA® Shoe cover, blue, embossed	41 x 15 cm	50/1000
211545	SELEFA® Shoe cover, blue, embossed	45 x 15 cm	50/1000

**Sales Unit:** Transport Carton

# Technical Summary

## Symbols



single use only



Contains natural rubber latex



### REACH regulation:

The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

### Compliance with legal requirements:

Medical Devices Directive 93/42/EEC, Class I

### Quality standards:

EN ISO 13485

### Label information and packaging compliance:

EN ISO 15233-1 and EN 1041

## Material

### Product part

### Raw Material & Characteristic

### Material

REF 211540 - CPE, white, 0.05 mm  
REF 211541 - CPE, blue, 0.035 mm  
REF 211545 - CPE, blue, 0.035 mm

## Warnings and Recommendations

### Warnings:

The product is single use. If reused performance of product may deteriorate, cross contamination may occur.

### Storage recommendation

Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.

### Disposal recommendation

Dispose in accordance with local regulations.


### Shelf life

3 Years

## Packing Information

Package	Material	Size ( L x W x H ) mm
Single pack	N/A	N/A
Dispenser pack	PE-bag	REF 211540 - 205 x 290 mm REF 211541 - 200 x 280 mm REF 211545 - 205 x 290 mm
Transport carton	Corrugated cardboard	REF 211540 - 400 x 295 x 390 mm REF 211541 - 390 x 282 x 368 mm REF 211545 - 390 x 290 x 375 mm
Packing tape	Transparent, PVC-free	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
211540	N/A	 6 438129 353037	 6 438129 553031
211541	N/A	 6 438129 352610	 6 438129 552614
211545	N/A	 6 438129 352740	 6 438129 552744