

# Gauze round sponges, non-sterile

REF: 222253, 222256

**SELEFA® Gauze round sponges, non-sterile** are highly absorbent sponges for a wide area of use such as: pre-operative skin cleansing, wound cleansing, drying, fluid absorption, protection etc.

## Key Features:

- ▶ 100% cotton gauze
- ▶ High absorption ability
- ▶ Skin and wound cleansing
- ▶ Fluid absorption



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
<b>222253</b>	Gauze round sponges, non-sterile	Ø 30	100/2000
<b>222256</b>	Gauze round sponges, non-sterile	Ø 70	100/1000

**Sales Unit:** Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices



**[ NATURAL LATEX FREE ]**

**REACH regulation:** The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-5 and EN ISO 10993-10

**Compliance with legal requirements:** Medical Devices Directive 93/42/EEC, Class I

**Compliance with product standards:** EN 14079

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1 and EN 1041

## Material

Product Part	Raw Material & Characteristic
Gauze	100% cotton gauze, 17 threads/cm <sup>2</sup> , white
Silicone ring	SEBS - Silicone Elastomer Blends, Latex free

## Warnings and Recommendations

**Warnings:** The product is single use. If reused performance of product may deteriorate, cross contamination may occur.

**Storage Recommendation** Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.

**Disposal Recommendation** Dispose in accordance with local regulations.

**Shelf Life** 5 Years

## Packing Information

Package	Material	Size ( L x W x H ) mm
Single Pack	N/A	N/A

<b>Dispenser Pack</b>	PE	222253: 250 x 350 222256: 400 x 550
<b>Transport Carton</b>	Carton, ECF or TCF quality	222253: 370 x 270 x 280 222256: 400 x 300 x 420
<b>Packing Tape</b>	Transparent, PVC-free	N/A

### Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>222253</b>	N/A	 6 438129 352290	 6 438129 552294
<b>222256</b>	N/A	 6 438129 352306	 6 438129 552300



Manufactured in China for  
**OneMed Group Oy**  
 Address: Metsäläntie 20  
 Postal Address: P.O. Box 10  
 FI-00321 Helsinki, Finland  
 Phone +358 20 786 6800  
 www.onemed.com  
 Domicile Helsinki  
 Business ID 2039640-1