

Gauze swabs, 2-pack, sterile

REF: 2222702, 2222722, 2222742, 2222762

SELEFA® Gauze swabs, 2-pack, sterile can be used for non-invasive medical procedures, such as wound padding, absorbing blood and exudates, skin and wound cleansing or as secondary dressing.

Key Features:

- ▶ Folded in cut edges to prevent loose threads
- ▶ 100% cotton gauze (17 threads per cm²)
- ▶ Sterilized by ethylene oxide gas (EO)
- ▶ Peroxide bleached (no chlorine)
- ▶ 2 swabs per pouch



REF	DESCRIPTION	SIZE	PACK COUNT unit/dispenser/carton
2222702	Gauze swabs, 2-pack, 8-ply, sterile	5 x 5 cm	2/120/2400
2222722	Gauze swabs, 2-pack, 8-ply, sterile	7,5 x 7,5 cm	2/120/2400
2222742	Gauze swabs, 2-pack, 8-ply, sterile	10 x 10 cm	2/120/2160
2222762	Gauze swabs, 2-pack, 12-ply, sterile	10 x 20 cm	2/40/720

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices
EO sterilized according to valid version of ISO 11135.



REACH regulation:

The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation:

EN ISO 10993-5, EN ISO 10993-7, EN ISO 10993-10

Compliance with legal requirements:

Medical Devices Directive 93/42/EEC, Class Is

Compliance with product standards:

EN 14079, EN ISO 11135-1, EN ISO 11737-1

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN 1041, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristic
--------------	-------------------------------

Main material	100 % cotton gauze, white
---------------	---------------------------

Warnings and Recommendations

Warnings:

The product is single use. If reused performance of product may deteriorate, cross contamination may occur.

Storage Recommendation

Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.

Disposal Recommendation

Dispose in accordance with local regulations.

Shelf Life

5 Years

Packing Information

Package	Material	Size (L x W x H) mm
---------	----------	-----------------------

Single Pack	Medical paper and Plastic film PP/PE	2222702: 72 x 88
		2222722: 98 x 113
		2222742: 122 x 138
		2222762: 133 x 249

Dispenser Pack	Carton, ECF or TCF quality	2222702: 90 x 70 x 180
		2222722: 115 x 95 x 180
		2222742: 140 x 120 x 180
		2222762: 250 x 130 x 80
Transport Carton	Carton, ECF or TCF quality	2222702: 370 x 190 x 370
		2222722: 495 x 240 x 370
		2222742: 435 x 370 x 370
		2222762: 510 x 400 x 250
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2222702	 6 438129 150971	 6 438129 350975	 6 438129 550979
2222722	 6 438129 150988	 6 438129 350982	 6 438129 550986
2222742	 6 438129 150995	 6 438129 350999	 6 438129 550993
2222762	 6 438129 151008	 6 438129 351002	 6 438129 551006