

SELEFA® Gauze swabs, 2-pack, 12-ply, sterile

REF: 222289, 222291

SELEFA® Gauze swabs, 12-ply, ES folded edges, sterile are primarily intended to cover or clean and absorb blood and exudate from wounds and abrasions. Can also be used for non-invasive medical procedures, such as wound padding, drying, protection, skin and wound cleansing or as secondary dressing. Good absorbent properties and in compliance with EN 14079 to assure high quality cotton gauze.

Key Features

- ▶ 100% cotton gauze
- ▶ Environmental friendly bleached (no chlorine)
- ▶ "ES" folded edges
- ▶ 17 threads per cm²
- ▶ Folded edges to prevent loose threads
- ▶ Sterilized using EO
- ▶ Easy-to-use unit pouch
- ▶ 2-pack




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
222289	Gauze swabs, 2-pack, 12-ply, sterile	10x10 cm	80 / 1280
222291	Gauze swabs, 2-pack, 12-ply, sterile, White	5x5 cm	80 / 1440

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class Is
Compliance with product standards:	EN 14079, EN ISO 11135-1, EN ISO 11737-1
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 20417, EN ISO 11607-1, EN ISO 15223-1

Material

Product Part	Raw Material & Characteristics
100 % cotton gauze, white	100 % cotton gauze, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	222289: 133 x 135 222291: 75 x 90
Dispenser Pack	TCF or ECF quality	222289: 135 x 130 x 145 222291: 90 x 78 x 145
Transport Carton	TCF or ECF quality	222289: 530 x 280 x 300 222291: 280 x 245 x 300
Packing Tape	PVC-free	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
222289	06438129151404	06438129351408	06438129551402
222291	6438129151428	6438129351422	6438129551426

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.