


SELEFA® Adhesive plaster, nonwoven

REF: 240100F, 240101F

SELEFA® Adhesive plaster, nonwoven is an injection plaster cut into strips. All with lined pad to prevent sticking to the wound.

Key Features

- Lined pad to prevent sticking to the wound




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
240100F	Adhesive plaster, nonwoven, Beige	2,4 x 7 cm	200 / 10200
240101F	Adhesive plaster, nonwoven, Beige	1,2 x 4 cm	200 / 10200

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-05 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class I
Compliance with product standards:	ASTM D3330
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristics
Carrier	Nonwoven, 100 % PES
Pad, surface	LDPE
Pad, absorbent core	Cellulose/PP/PES
Skin adhesive	Synthetic, colofonium-free
Release paper	Silicone

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Carton, TCF or ECF quality	93 x 165 x 35
Transport Carton	Carton, TCF or ECF quality	240100F: 590 x 390 x 430 240101F: 580 x 185 x 300
Packing Tape	PVC-free, Transparent	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
240100F	N/A	6438129352177	6438129552171
240101F	N/A	6438129352184	6438129552188

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.