

evercare[®] Optimia SPECIAL-OP-gown

REF: 2576-01, 2577-01, 2578-01, 2579-01, 2580-01 with hand towels

evercare[®] Optimia SPECIAL-OP-gown is a SMS gown designed to provide ultimate comfort and protection. The gown aim to provide a cost efficient alternative to the traditional gowns made of cotton like fabrics. The gown complies to the EN13795 for surgical gowns, high performance.

Key Features:

- ▶ Soft SMS material with alcohol repellent treatment
- ▶ Breathable reinforcement
- ▶ Ergonomic design in shoulders for freedom of movement
- ▶ Long arms to provide comfort and reduce tensions
- ▶ Belts positioned to improve comfort
- ▶ The adjustment feature in neck is sewed in a double stitching
- ▶ Suitable for invasive surgery



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2576-01	evercare [®] Optimia SPECIAL-OP-gown, with hand towels	M-short, length 120 cm	20/40
2577-01	evercare [®] Optimia SPECIAL-OP-gown, with hand towels	M, length 130 cm	18/36
2578-01	evercare [®] Optimia SPECIAL-OP-gown, with hand towels	L, length 140 cm	18/36
2579-01	evercare [®] Optimia SPECIAL-OP-gown, with hand towels	XL, length 150 cm	16/32
2580-01	evercare [®] Optimia SPECIAL-OP-gown, with hand towels	XXL, length 160 cm	16/32

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices
EO sterilized according to valid version of ISO 11135



REACH regulation:	The product covered by this product data sheet does not contain carcinogenic, mutagenic or toxic substances, including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	ISO 10993-5, ISO 10993-10
Compliance with legal requirements:	Medical Devices Directive 93/42/EEC, Class Is
Compliance with product standards:	EN 13795
Quality standards:	ISO13485
Label information and packaging compliance:	EN ISO 15223-1 and EN 1041

Material

Product Part	Raw Material & Characteristic
Main material	SMS nonwoven, multilayer, PP, blue
Reinforcement sleeves and front	Breathable laminate made of PP nonwoven and PE-film
Cuffs	Rib knitted fabric, PES, white
Neck ribbon	PP nonwoven, white
Velcro band	Woven ribbon with hook/loop, PA
Belt card	with evercare® printing, pulp(ECF), white
Sewing thread	PES, white
Hand towel	Natural white paper thread reinforced, pulp(ECF), PES(reinf.threads)
Crepe wrapping	Natural white paper, pulp(ECF)

Warnings and Recommendations

Warnings: The product is single use. If reused performance of product may deteriorate, cross contamination may occur.

Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight , other intensive light sources and ozone.
Disposal Recommendation	Product and packaging materials to energy recovery or mixed waste.
Shelf Life	5 Years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Paper(TCF)	2576-01: 250 x 390
		2577-01: 250 x 390
		2578-01: 250 x 390
		2579-01: 280 x 415
		2580-01: 280 x 415
Dispenser Pack	Carton, TCF quality	381 x 279 x 320
Transport Carton	Carton, TCF quality	574 x 401 x 330
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2576-01	 6 438129 123586	 6 438129 323580	 6 438129 523584
2577-01	 6 438129 123593	 6 438129 323597	 6 438129 523591
2578-01	 6 438129 123609	 6 438129 323603	 6 438129 523607
2579-01	 6 438129 123616	 6 438129 323610	 6 438129 523614
2580-01	 6 438129 123623	 6 438129 323627	 6 438129 523621



OneMed Group Oy
 Address: Metsäläntie 20
 Postal Address: P.O. Box 10
 FI-00321 Helsinki, Finland
 Phone +358 20 786 6800
 www.onemed.com
 Domicile Helsinki
 Business ID 2039640-1

