

evercare® Round cap Rondo Standard

REF: 2740-02, 2741-02, 2742-02, 2750-02, 2751-02, 2752-02, 2760-02, 2761-02, 2762-02

evercare® round cap Rondo Standard is made of a new light weight SMS material that provides high protection against particles from the scalp of the wearer as well as good fluid protection. The lower use of polypropylene material means a reduced environmental impact. Rondo is an economical range of headwear covering all basic needs. Rondo is also worn by medical personnel and patients in other environments for hygiene reasons. Available in different colors and sizes. The use of headwear is proven to reduce the risk of particles falling into the surgical wound. * Fits the wall dispenser REF 2900

Key Features:

- ▶ Reduces the risk of particles falling into the wound
- ▶ Reduced environmental impact
- ▶ M, L and XL size
- ▶ Different colors
- ▶ Fits the wall dispenser REF 2900




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2740-02	Round cap Rondo Standard, green	M (Ø 48 cm)	150/900
2741-02	Round cap Rondo Standard, white	M (Ø 48 cm)	150/900
2742-02	Round cap Rondo Standard, blue	M (Ø 48 cm)	150/900
2750-02	Round cap Rondo Standard, green	L (Ø 52 cm)	150/900
2751-02	Round cap Rondo Standard, white	L (Ø 52 cm)	150/900
2752-02	Round cap Rondo Standard, blue	L (Ø 52 cm)	150/900
2760-02	Round cap Rondo Standard, green	XL (Ø 55 cm)	150/900
2761-02	Round cap Rondo Standard, white	XL (Ø 55 cm)	150/900
2762-02	Round cap Rondo Standard, blue	XL (Ø 55 cm)	150/900

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices

**REACH regulation:**

The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation:

EN ISO 10993-05, EN ISO 10993-10

Compliance with legal requirements:

Medical Devices Directive 93/42/EEC, Class I

Compliance with product standards:

N/A

Quality standards:

EN ISO 13485

Label information and packaging compliance:

EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristic
Main material	PP SMS nonwoven
Elastic band	Yarn covered synthetic rubber, natural latex free
Sewing thread	Polyester, white

Warnings and Recommendations

Warnings:

The product is single use. If reused performance of product may deteriorate, cross contamination may occur.

Storage Recommendation

Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.

Disposal Recommendation

Dispose in accordance with local regulations.



















Shelf Life

5 Years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	Carton, TCF or ECF quality	280 x 270 x 180
Transport Carton	Carton, TCF or ECF quality	560 x 295 x 555
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2740-02	N/A	 6 438129 320756	 6 438129 520750
2741-02	N/A	 6 438129 317978	 6 438129 520743
2742-02	N/A	 6 438129 320756	 6 438129 520750
2750-02	N/A	 6 438129 320787	 6 438129 520781
2751-02	N/A	 6 438129 318036	 6 438129 520798
2752-02	N/A	 6 438129 320800	 6 438129 520804
2760-02	N/A	 6 438129 320831	 6 438129 520835
2761-02	N/A	 6 438129 320848	 6 438129 520842
2762-02	N/A	 6 438129 320855	 6 438129 520859