

evercare® Round cap Rondo Standard

REF: 2740-02B, 2742-02, 2750-02, 2751-02, 2752-02, 2760-02, 2761-02, 2762-02

evercare® Headweares are designed to help prevent foreign particles such as hair or skin scales falling from medical personnel into the wound during surgical interventions. evercare® round cap Rondo Standard is made of a new light weight SMS material that provides high protection against particles from the scalp of the wearer as well as good fluid protection. The lower use of polypropylene material means a reduced environmental impact. Rondo is an economical range of headwear covering all basic needs. Rondo is also worn by medical personnel and patients in other environments for hygiene reasons. Available in different colors and sizes.

Key Features

- ▶ Reduces the risk of particles falling into the wound
- ▶ Reduced environmental impact
- ▶ M, L and XL size
- ▶ Different colors




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2740-02B	Round cap Rondo Standard, Green	M, 48 cm	150 / 900
2742-02	Round cap Rondo Standard, Blue	M, 48 cm	150 / 900
2750-02	Round cap Rondo Standard, Green	L, 52 cm	150 / 900
2751-02	Round cap Rondo Standard, White	L, 52 cm	150 / 900
2752-02	Round cap Rondo Standard, Blue	L, 52 cm	150 / 900
2760-02	Round cap Rondo Standard, Green	XL, 55 cm	150 / 900
2761-02	Round cap Rondo Standard, White	XL, 55 cm	150 / 900
2762-02	Round cap Rondo Standard, Blue	XL, 55 cm	150 / 900

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-05 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDD class I
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristics
Main material	PP SMS nonwoven
Elastic band	Yarn covered synthetic round rubber, natural latex free
Sewing thread	Polyester, white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
---------	----------	-----------------------

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2740-02B	N/A	06438129320732	06438129520736
2742-02	N/A	06438129320756	06438129520750
2750-02	N/A	6438129320787	6438129520781
2751-02	N/A	6438129320794	6438129520798
2752-02	N/A	06438129320800	06438129520804
2760-02	N/A	6438129320831	6438129520835
2761-02	N/A	6438129320848	6438129520842
2762-02	N/A	06438129320855	06438129520859

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.