



## evercare® Surgical Hood William

REF: 2786

evercare® Headweares are designed to help prevent foreign particles such as hair or skin scales falling from medical personnel into the wound during surgical interventions. These very popular surgical hoods come with different special features and designs, such as integrated sweatbands for additional comfort, or open design for easier dressing placed on the head instead of over the head. evercare® Surgical Hood William is a new low lint hood especially designed for surgical staff. The new tape adjustment feature in neck makes it easier to wear a headlamp. Integrated sweatband, shoulder protection and soft binding around the cheek is included and provides extra comfort. The headwear material complies to the EN 13795 for clean air suit to contribute to ultra clean air. The product is separately packed to provide extra cleanliness. Users who get skin irritation and are sensitive for viscose material have switched to William.

### Key Features

- ▶ Reduces the risk of particles falling into the wound
- ▶ Universal size
- ▶ Shoulder protection
- ▶ Integrated sweat band
- ▶ Complies to the EN 13795 for clean air suit
- ▶ The product is separately packed to provide extra cleanliness




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2786	Surgical hood William, shoulder protection, sweatband, Blue	Universal	30 / 180

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Compliance with product standards:</b>	EN 13795-2 for Clean Air Suit
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
<b>Main material</b>	PP SMS nonwoven, blue
<b>Sweat band, face lining</b>	Viscose spunlance nonwoven, white
<b>Elastic band</b>	Nylon and synthetic rubber, latex free, white
<b>Sewing thread</b>	Nylon and spun polyester, white
<b>Fixation tapes</b>	Double adhesive tape, latex-free, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	2,5 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	PE bag, Transparent	
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	300 x 185 x 145
<b>Transport Carton</b>	Carton, TCF or ECF quality	385 x 310 x 455
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>2786</b>	N/A	6438129319569	6438129519563

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.