

# evercare® Medical face mask, with visor and tie-bands, type IIR

REF: 2854

evercare® Medical face mask, Face shield with tie-bands is a type IIR face mask with visor that provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is a splash resistant mask offering protection against potentially contaminated fluid splashes. Recommended to be used during wet surgery. Appropriate for situations where exposure to blood and/or bodily fluids from the patient or work environment is a risk for the healthcare worker.

## Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Splash resistant
- ▶ Color: Outer facing blue, inner facing white
- ▶ Anti-fog for a clear view
- ▶ Soft inner layer for comfort
- ▶ Anti-reflection to prevent glare



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2854	Medical face mask, with anti-reflective visor, anti-fog, tie-bands, type IIR, non-sterile., Blue	178 x 95 mm	25 / 100

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-05 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Compliance with product standards:</b>	EN 14683:2019, Type IIR, Bacterial Filtration Efficiency (BFE) ≥98%
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
<b>1st outer coverstock</b>	Spunbond PP, blue
<b>2nd outer coverstock</b>	Spunbond PP, white
<b>Filter</b>	Meltblown PP, white
<b>Inner coverstock</b>	Hydrophobic tissue, white
<b>Upper binding</b>	Rayon and PET drylaid, black
<b>Lower binding</b>	Spunbond PP, white
<b>Ties</b>	100 % PES chemicalbonded, white, upper/lower length 433 mm
<b>Nose wire</b>	Pe-coated galvanized wire, 125 mm
<b>Plastic shield</b>	Anti-fog PET film, transparent
<b>Foam</b>	PU foam, blue + Spunbond PP, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	3 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Corrugate Paper, Flute E	220 x 110 x 130
<b>Transport Carton</b>	Carton, TCF or ECF quality	345 x 232 x 280
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>2854</b>	N/A	6438129315639	6438129515633

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.