

# evercare<sup>®</sup> Medical face mask Multiguard Antifog, special, with tie-bands, type II

REF: 2865

evercare<sup>®</sup> Medical face mask Multiguard Antifog, special with tie-bands type II provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is an allround face mask with high filtration capacity and low breathing resistance. Recommended to be used during dry surgery, with no risk of splashes.

## Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Color: Outer facing green, inner facing white
- ▶ Antifog



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2865	Medical face mask Multiguard Antifog Special with tie-bands, type II, non-sterile, Green	175 x 90 mm	50 / 300

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-05 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Compliance with product standards:</b>	EN 14683:2019, Type II, Bacterial Filtration Efficiency (BFE) ≥98%
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
<b>Outer facing</b>	Spunbond PP, green
<b>Filter media</b>	Meltblown PP, white
<b>Inner facing</b>	Thermalbonded PP, white
<b>Bindings</b>	Spunbond PP, white, upper/lower length 435 mm
<b>Nose wire</b>	PE coated galvanize wire, 125 mm
<b>Anti-fog front</b>	Spunbond PP, white + Coated LDPE, opaque white
<b>Anti-fog back</b>	PU foam, light blue + Spunbond PP, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	3 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
Single Pack	N/A	N/A
Dispenser Pack	Coated paper	220 x 110 x 130
Transport Carton	Carton, TCF or ECF quality	345 x 232 x 280
Packing Tape	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
2865	N/A	6438129312973	6438129512977

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.