

evercare® Medical face mask with earloops, type II

REF: 2882,2883

evercare® Medical face mask with ear-loops, type II provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. The mask is slightly shorter in width than many other face masks but at the same time it is also wider when in open size. The design of the side sealing and the anchoring of the ear loops help to reduce the gaping at the cheeks and makes the mask more conforming to the face. This overall design results in a deeper cup-shape, which is better conforming and provides a fitter coverage of the nose and mouth. This is a single use product.

Key Features:

- ▶ Comfortable elastic earloops design for reduced pressure to the ears
- ▶ Anatomic pleated folding limiting inward leakage
- ▶ EN 14683 type II 3ply face mask with ear-loops
- ▶ Designed to be more conforming to the face
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Flexible nose wire for proper fit
- ▶ Color: White



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2882	Medical face mask with earloops, type II, white	165 x 90 mm	50/2500
2883	Medical face mask with earloops, type II, 10-pack, white	165 x 90 mm	120/2160

Sales Unit: Transport Carton

Technical Summary

Regulations

The product complies with legal requirements for medical devices.



REACH regulation:	The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-05, EN ISO 10993-10
Compliance with legal requirements:	Medical Devices Directive 93/42/EEC, Class I
Compliance with product standards:	EN 14683:2019, Type II, Bacterial Filtration Efficiency (BFE) ≥98%
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

Material

Product Part	Raw Material & Characteristic
Outer facing	Spunbond PP nonwoven, white
Filter media	Meltblown PP nonwoven, white
Inner facing	Spunbond PP nonwoven, white
Elastic earloops	Spandex, Polyester & Polyamide, latex-free, white
Nose wire	100% PE

Warnings and Recommendations

Warnings:	The product is single use. If reused performance of product may deteriorate, cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	2882: N/A 2883: PE-bag, transparent	2882: N/A 2883: 125 x 212
Dispenser Pack	Carton, TCF or ECF quality	2882: 188 x 98 x 87 2883: 190 x 123 x 188
Transport Carton	Carton, TCF or ECF quality	2882: 515 x 395 x 460 2883: 590 x 390 x 400
Packing Tape	Transparent, PVC-free	N/A

Barcodes

REF	Single Pack	Dispenser	Transport Carton
2882	N/A	 6 438129 324624	 6 438129 524628
2883	 6 438129 124637	 6 438129 324631	 6 438129 524635



OneMed Group Oy
 Address: Metsäläntie 20
 Postal Address: P.O. Box 10
 FI-00321 Helsinki, Finland
 Phone +358 20 786 6800
 www.onemed.com
 Domicile Helsinki
 Business ID 2039640-1