


evercare[®] Medical Face Mask typellR, tieband, antifog, eye shield, blue

REF: 2886

evercare[®] Medical face mask, Face shield with tie-bands is a type IIR face mask with visor that provides a barrier to minimize the direct transmission of infective agents between medical staff and patient. It is a splash resistant mask offering protection against potentially contaminated fluid splashes. Recommended to be used during wet surgery. Appropriate for situations where exposure to blood and/or bodily fluids from the patient or work environment is a risk for the healthcare worker.

Key Features

- ▶ EN 14683:2019
- ▶ Bacterial Filtration Efficiency (BFE) ≥ 98%
- ▶ Single use
- ▶ Flexible nose wire for proper fit
- ▶ Splash resistant
- ▶ Anti-fog for a clear view
- ▶ Soft inner layer for comfort
- ▶ Anti-reflection to prevent glare
- ▶ Color: Outer facing blue, inner facing white




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2886	evercare [®] Medical Face Mask typellR, tie-band, antifog, eye shield, blue, Blue	17,5 x 9,5 cm	25 / 100

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-5
EN ISO 10993-10

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR class I

Compliance with product standards: EN14683:2019

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 20417

Material

Product Part	Raw Material & Characteristics
1st Outer facing	Spunbond PP, blue
2nd Outer facing	Spunbond PP, white
Filter media	Meltblown PP, white
Inner facing	Hydrophobic tissue, white
Tie-bands	PP Spoundbond, white
Nose wire	PE coated metal wire, 125 mm
Plastic shield	Anti-fog PET film, transparent
Foam	PU foam, grey

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 Years

Packing information

Package	Material	Size (L x W x H) mm
Single Pack	N/A	N/A
Dispenser Pack	TCF or ECF quality	222 x 190 x 90
Transport Carton	TCF or ECF quality	392 x 323 x 200
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
2886	N/A	06438129325089	06438129525083

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.