

## evercare<sup>®</sup> Surgical swabs, X-ray detectable, nonwoven, sterile







REF: 3230, 3231, 3232, 3233, 3234, 3235

evercare<sup>®</sup> Surgical swabs X-ray detectable, nonwoven, sterile have color coding for easy product identification and numeric control system for counting assistance. The control system contains following information: color code, REF number, dry weight of the product, number of individual swabs, size of individual swabs, material, and control batch number. As part of the production process the products are counted and controlled by a metal detector and the material is vacuumed to minimize loose threads. These swabs are low linting and cost efficient. Includes extra nonwoven pouch under the sterile barrier, to allow aseptic technique when taking out from the pouch and keeping the products in place.

### Key Features

- ▶ X-ray detectable
- ▶ Sterilized with EO
- ▶ Numeric control system
- ▶ Colour coding system for easy product identification
- ▶ Wrapped into a sealed lint-free easy-to-open wrap
- ▶ High absorbency
- ▶ X-ray detectable thread safely welded to the nonwoven




REF	DESCRIPTION	COLOR CODING	SIZE	Pack Count dispenser/carton
3230	Surgical swabs, X-ray detectable, nonwoven, 5-pack, sterile, Orange		10x40 cm	130 / 1040
3231	Surgical swabs, X-ray detectable, nonwoven, 5-pack, sterile, Green		20x40 cm	65 / 520
3232	Surgical swabs, X-ray detectable, nonwoven, 5-pack, sterile, Brown		30x40 cm	45 / 360
3233	Surgical swabs, X-ray detectable, nonwoven, single pack, sterile, Grey		40x65 cm	20 / 160
3234	Surgical swabs, X-ray detectable, nonwoven, 5-pack, sterile, Black		40x65 cm	20 / 160
3235	Surgical swabs, X-ray detectable, nonwoven, 5-pack, sterile, Blue		40x40 cm	30 / 240

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



**REACH regulation:** The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-5  
EN ISO 10993-7  
EN ISO 10993-10  
EN ISO 10993-11

**Compliance with legal requirements:** Medical Device Regulation (EU) 2017/745, MDR class IIa

**Compliance with product standards:** EN ISO 11135, EN 1644-1, EN 1644-2

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
Nonwoven	Viscose 70 % Polyester 30 %, white
X-ray detectable thread	DEHP-free, black X-ray thread, with min. 60 % (w/w) barium sulphate
Inner pouch	PP nonwoven, blue

## Warnings & Recommendations

**Warnings** The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.  
Control the amount of items before use. If the amount is different than stated on the packaging, do not use the product.

**Storage Recommendation** Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.

**Disposal Recommendation** Dispose in accordance with local regulations.

**Shelf Life** 5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, PP/PE film	3230: 200 x 185
		3231: 200 x 185
		3232: 200 x 185
		3233: 200 x 265
		3234: 200 x 165
		3235: 200 x 155
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	3230: 190 x 165 x 285
		3231: 190 x 165 x 285
		3232: 190 x 165 x 285
		3233: 285 x 175 x 190
		3234: 285 x 175 x 190
		3235: 190 x 165 x 285
<b>Transport Carton</b>	Carton, TCF or ECF quality	3230: 600 x 400 x 350
		3231: 600 x 400 x 350
		3232: 600 x 400 x 350
		3233: 600 x 400 x 370
		3234: 600 x 400 x 370
		3235: 600 x 400 x 350
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
3230	06438129114362	06438129314366	06438129514360
3231	06438129114379	06438129314373	06438129514377
3232	06438129114386	06438129314380	06438129514384
3233	06438129114393	06438129314397	06438129514391
3234	06438129114409	06438129314403	06438129514407
3235	06438129119022	06438129319026	06438129519020

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.