


evercare[®] Medikit Washing set 590B

REF: 400607


evercare[®] MediKit Single-use procedure small kits are assembled to provide healthcare professionals with the necessary supplies for specific medical interventions. The primary purpose of the washing kit is to facilitate wound management tasks effectively, when it involves cleansing and irrigating of wounds. Available in a single, sterile package, ensuring optimal hygiene conditions and patient safety.

Key Features

- ▶ 30 Nonwoven swabs 5x5cm 30g 4ply
- ▶ 1 Nonwoven swabs 7,5x7,5cm 30g 4ply
- ▶ Pouch bottom side has tray function



30 Nonwoven swabs
5 x 5 cm, 30 g, 4 ply



1 Nonwoven swab
7,5 x 7,5 cm, 30 g, 4 ply

REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
400607	Washing Set 590B		66 / 132

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-7

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR Class Is

Compliance with product standards: EN ISO 11135

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 20417, EN ISO 11607-1, EN ISO 15223-1

Material

Product Part	Raw Material & Characteristics
30 Nonwoven swabs 5x5cm 30g 4ply	VIS+PES white
1 Nonwoven swabs 7,5x7,5cm 30g 4ply	VIS+PES white

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packaging Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	100 x 200 x 30
Dispenser Pack	Cardboard	380 x 288 x 335
Transport Carton	TCF or ECF quality	595 x 400 x 355
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser Pack	Transport Carton
400607	06942013200890	16942013200897	26942013200894

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.