

470

Clip fasteners

Grampos para ligaduras

Agrafes pour pansement

Verbandklammern

Clipsluiting

Graffette

Grapas para vendas

Clips



Clip fasteners

- With metallic ends for bandage fixation.
- Easily repositionated.
- Strong fixation.
- Can be used on all types of bandages.
- Shelf box + transport carton.

Non sterile

Code	Color	Pcs./Sh.box	Pcs./Tr. carton	EAN Sh.box	EAN Tr. carton
470-001	White/Branco	100	6.400	560562240831	560562200831
470-002	Cream/Creme	100	6.400	560562240832	560562200832