

# PRODUCT SPECIFICATION



Art.No: 500121 Protective sheet 37x50

<b>End use:</b>	Protection against liquids.																														
<b>Description:</b>	<p>The sheets is a laminate made of one layer of highly absorbent tissue and PE-film (polyethylene) 37x50 cm.</p> <p>The products are packed in plastic bags and then in cartons on pallets. The pallets are stretch wrapped.</p>																														
<b>Material, product:</b>	<p>The tissue is made of oxygen bleached cellulose fibres, TCF-pulp (total chlorine free), 100 % virgin fibres.</p> <p>The film is blue LDPE (polyethylene).</p> <table><tr><td>Tissue basis weight:</td><td>21,0</td><td>gsm</td></tr><tr><td>PE-film: Blue</td><td>11</td><td>my</td></tr></table>	Tissue basis weight:	21,0	gsm	PE-film: Blue	11	my																								
Tissue basis weight:	21,0	gsm																													
PE-film: Blue	11	my																													
<b>Packaging:</b>	<table><tr><td>Units/inner package:</td><td>200</td><td>pcs</td></tr><tr><td>Inner package/transport units:</td><td>6</td><td>pcs</td></tr><tr><td>Transport units/pallet:</td><td>30</td><td>pcs</td></tr><tr><td>Weight, transport unit</td><td>Net</td><td>7,1 kg</td></tr><tr><td></td><td>Gross:</td><td>7,5 kg</td></tr><tr><td>Measurement/transport unit</td><td>Length:</td><td>390 mm</td></tr><tr><td></td><td>Width:</td><td>390 mm</td></tr><tr><td></td><td>Height:</td><td>245 mm</td></tr><tr><td>Pallet:</td><td>EUR</td><td>800x1200 mm</td></tr><tr><td></td><td>Height:</td><td>1375 mm</td></tr></table>	Units/inner package:	200	pcs	Inner package/transport units:	6	pcs	Transport units/pallet:	30	pcs	Weight, transport unit	Net	7,1 kg		Gross:	7,5 kg	Measurement/transport unit	Length:	390 mm		Width:	390 mm		Height:	245 mm	Pallet:	EUR	800x1200 mm		Height:	1375 mm
Units/inner package:	200	pcs																													
Inner package/transport units:	6	pcs																													
Transport units/pallet:	30	pcs																													
Weight, transport unit	Net	7,1 kg																													
	Gross:	7,5 kg																													
Measurement/transport unit	Length:	390 mm																													
	Width:	390 mm																													
	Height:	245 mm																													
Pallet:	EUR	800x1200 mm																													
	Height:	1375 mm																													
<b>Material, packaging:</b>	<p>The bags are made of transp PE (polyethylene) produced of rawmaterial that fulfills the requirements for food contacts as per BGA and FDA.</p> <p>The cartons are made of brown cardboard.</p> <p>The stretch film is made of transparent PE (polyethylene) film.</p>																														
<b>Destruction:</b>	<p>The cartons could be incinerated and recycled as cardboard.</p> <p>The inner packaging and stretch film could be incinerated under standard conditions or recycled as polyethylene film (soft plastic)</p>																														
<b>Certification:</b>	<p>Finess Hygiene AB is certified in accordance with ISO 9001:2008 13485:2012 and ISO 14001:2004. Our CE-approved products fulfil the requirements of Directive MDD 93/42/EEC.</p>																														

Kisa 161028

**Finess Hygiene AB** SE-590 40 KISA  
Tel:+46(0)49415100 Fax:+46(0)49415191