

bedal pivc

Suitable for stabilizing & fixating a large range of catheters with a luer lock connection, such as: Peripheral IV lines, arterial lines & mini-midlines.

✓ **Reducing risk of catheter dislodgement (migration, pull out)**

- High fixation strength
-

✓ **Reducing infection risk**

- Absence of catheter-pistoning (micro motions)
 - Absence of any glue to fixate catheter onto device, which reduces dirt build-up
 - Tape material is non-absorbent and resistant to excess fluids
-

✓ **Easy & time-saving replacements**

- Easy & safe separating of PU dressing from **bedal pivc**, due to 'non-stick' plastics & textured surface of tape
-

✓ **Easy catheter cleaning**


- Absence of any glue to fixate catheter onto device; no glue residue on catheter
-

✓ **Optimal patient comfort**

- Tape is skin friendly (hypoallergenic), highly breathable with foam tape providing cushioning effect on skin
 - Low profile
 - Use of soft materials only
-

✓ **Easy to use**

- Simple functioning of the device; opening & closing of the leash
 - Catheter re-positioning possible by re-opening & closing the leash
 - Easy handling & placement of device while wearing gloves
-



Ordering information

Product number	Description	Size range	Qty
PV02	bedal pivc	all catheters with luer lock	25 pcs / box

Find out more on www.bedal.be