

TP-MF-B



IV Accessories

SKU/REF TP-MF-B

SKU/REF Name Male-Female Luer Lock Cap, white

Resources (1)

 [Infusion specialty disposables e-catalogue](#)

Specifications

Country of Manufacture	China
Classification	Is
Quantity per Box	100
Quantity per Case	1000
Colour Code	White
Shelf Life	48 Months
Sterilization Method	Ethylene Oxide
PVC	PVC-Free
Material	Polyethylene (PE)
DEHP	DEHP-Free

Please note, not all products, services or features of products and services may be available in your local area. Please check with your local BD representative.



COMPANY

[About BD](#)

[BD values](#)

[Regulatory documents](#)

[Careers](#)

SUPPORT

[Support](#)

[Contact us](#)

[COVID-19](#)



© 2023 BD. All rights reserved. Unless otherwise noted, BD, the BD Logo and all other trademarks are property of Becton, Dickinson and Company.

[Privacy Policy](#) | [Terms of Use](#) | [Cookie-inställningar](#)