

bedal fixation devices

The Bedal portfolio of fixation devices: accelerating in versatility:

bedal flex	Plexus & epidural catheters, spinal neuro modulation electrodes, harpoons, CVC & midlines.
bedal air	Drainage catheters, CVC, Gripper needles, thrombolysis catheters, tubings, cables & drivelines.
bedal picc cvc	All types & brands of venous catheters.
bedal pivc	Peripheral IV lines, arterial lines & mini-midlines

✓ Reducing risk of catheter dislodgement (migration, pull out)

- High fixation strength

✓ Reducing infection risk

- Absence of catheter-pistoning (micro motions)
- Anti-leakage snap-part prevents body fluids & dirt build-up to enter into device
- Absence of any glue to fixate catheter onto device, which reduces dirt build-up
- Tape material is non-absorbent and resistant to excess fluids

✓ Easy & time-saving replacements

- Easy & safe separating of PU dressing from **bedal** products, due to 'non-stick' plastics & textured surface of tape

✓ Easy catheter cleaning



- Absence of any glue to fixate catheter onto device; no glue residue on catheter

✓ Optimal patient comfort

- Tape is skin friendly (hypoallergenic), highly breathable with foam tape providing cushioning effect on skin
- Use of soft materials

✓ Easy to use

- Simple functioning of the device; opening & closing of the leash
- Catheter re-positioning possible by re-opening & closing the leash
- Easy handling & placement of device while wearing gloves



Ordering information

Product number	Description	Size range	Qty
EP02	bedal flex , size XS	2.7-5 Fr / 20-16 G	50 pcs / box
UN02	bedal air , size S	3-12 Fr	25 pcs / box
DP01	bedal air , size M	6-16 Fr	25 pcs / box
PC02	bedal picc cvc	all picc/cvc catheters	25 pcs / box
PV02	bedal pivc	all catheters with luer lock	25 pcs / box

Find out more on www.bedal.be