

# Ryles Tube Products (ENFit)

## Product Data Sheet

<b>Description</b>	Ryles Tube, Graduated, X-Ray with closed tip, four side eyes and ENFit connector.
<b>Intended Use</b>	Ryles Tube is designed to provide access to nasogastric tract for introduction of nutrition and aspiration of intestinal secretion. The main use of a Nasogastric tube is for feeding and for administering drugs and other oral agents. A Nasogastric tube is inserted through the nose with the end of the tube reaching into the stomach.
<b>Duration of Use</b>	Short Term Use: for more than 60 minutes and less than 30 days.
<b>Materials</b>	Tube - PVC (DEHP Free); ENFit Connector - ABS; Loop - LDPE; Radio-Opaque Line (X-Ray). Not made with DEHP Plasticizer. Not made with Natural Rubber Latex.
<b>Product Information</b>	Siliconised and smooth outer surface to reduce chances of tissue irritation. Soft, frosted and kink resistance PVC tubing. 105cm tube with radio-opaque line provided to its full length and marked at 50, 60, and 70 cm from the tip for accurate placement of tube. Closed distal end with four lateral eyes to reduce risk of blockage. ENFit Enteral Feeding Connector - has a unique enteral specific design in accordance to ISO 80369-3 that does not allow connectivity with connectors for any other clinical use. It provides a locking feature for a secure, leak free connection. Tube shore hardness: 75~78A. Available in different sizes from FG 6 to FG 18. Supplied sterile, in individual peel pouches.
<b>Quality System</b>	EN ISO 13485:2016
<b>EU Classification</b>	Class IIa
<b>Sterilised by</b>	Ethylene Oxide
<b>Shelf Life</b>	5 years
<b>Storage &amp; Safety</b>	Toxicity: Manufactured from biocompatible materials. Storage: Store in between 5°C to 35°C, avoid excessive heat, protect from direct sunlight and moisture.


Product Code	Tube Length	CH/Dia	Connector	Side Eyes	Tip	Box Qty
AERT2306EF	105cm	6ch ø 2.0mm	ENFit	4	Closed Tip	50
AERT2308EF	105cm	8ch ø 2.7mm	ENFit	4	Closed Tip	50
AERT2310EF	105cm	10ch ø 3.3mm	ENFit	4	Closed Tip	50
AERT2312EF	105cm	12ch ø 4.0mm	ENFit	4	Closed Tip	50
AERT2314EF	105cm	14ch ø 4.7mm	ENFit	4	Closed Tip	50
AERT2316EF	105cm	16ch ø 5.3mm	ENFit	4	Closed Tip	50
AERT2318EF	105cm	18ch ø 6.0mm	ENFit	4	Closed Tip	40

**\*Please note Pennine does not recommend non ENFit tubes for feeding in markets that have adopted the ENFit System.**

© 2022 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.


 Poly Medisure Limited, Plot No. 104-105, Sector-59, HSIIDC Industrial Area, Ballabhgarh, Faridabad, Haryana 121004, INDIA

**EC REP** Obelis S.A. Boulevard General Wahis 53, B-1030, Brussels, Belgium. mail@obelis.net



# Ryles Tube Products (ENFit)

## Product Data Sheet




\*Please note Pennine does not recommend non ENFit tubes for feeding in markets that have adopted the ENFit System.

© 2023 Pennine Healthcare or one of its subsidiaries.  
All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

 Poly Medisure Limited, Plot No. 104-105, Sector-59, HSIIDC Industrial Area, Ballabhgarh, Faridabad, Haryana 121004, INDIA

**EC REP** Obelis S.A. Boulevard General Wahis 53, B-1030, Brussels, Belgium. mail@obelis.net

**STERILE EO MD**  **CE** **0123**   