


Link Yankauers

Product Data Sheet

Product Codes	See table over page
Description	Single use, rigid suction devices for use with suction connecting tubing for general purpose surgical suctioning.
Intended Use	Surgical suction probe range of devices are hand operated, single use, and are used as a conduit to remove bodily fluids and secretions, surgical tissue debris, irrigation fluids and gasses. They are used in surgical sites and body cavities to facilitate observation and/or to clear an obstruction. To be used on all patient populations by clinicians and trained homecare users, in clinical setting, homecare and emergency settings.
Duration of Use	Normally intended for continuous use of less than 60 minutes, but will be patient specific depending on the duration of the surgical procedure.
Product Information	Available with or without vacuum control (VC) vent Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip
Materials	PVC (DEHP Free) Latex Free
Quality System	ISO 13485
EU Device Classification	Class IIa
Box Qty	See table over page
Packaging	Sterile, double wrapped individual peel pouches
Sterilised by	Ethylene Oxide;
Shelf Life	5 years
Storage & Safety	Toxicity: Manufactured from biocompatible materials. Storage: Store in a clean, cool, dry place away from direct sunlight



© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

Link Yankauers

Product Data Sheet

Product Code	Type	Length	CH	Qty per box
LY-3400	Short Mini Tip	190mm	15ch	100
LY-3401	Short Mini Tip VC	190mm	15ch	100
LY-4400	Short Mini Super Tip	195mm	15ch	100
LY-4401	Short Mini Super Tip (VC)	195mm	15ch	100
LY-4420	Short Mini Super Tip (2 side)	195mm	15ch	100
LY-3500	Mini Fine Tip	255mm	12ch	50
LY-3501	Mini Fine Tip (VC)	255mm	12ch	50
LY-4500	Mini Fine Super Tip	255mm	12ch	50
LY-4500/NL	Mini Fine Super Tip (INT)	255mm	12ch	50
LY-4501	Mini Fine Super Tip (VC)	255mm	12ch	50
LY-4501/NL	Mini Fine Super Tip (VC) (INT)	255mm	12ch	50
LY-3600	Midi Tip	240mm	18ch	50
LY-3601	Midi Tip VC	240mm	18ch	50
LY-3600/PR	Midi Tip (PR)	240mm	18ch	50
LY-4600	Midi Super Tip	240mm	18ch	50
LY-4600/PR	Midi Super Tip (PR)	240mm <td 18ch	50	
LY-4601	Midi Super Tip (VC)	240mm	18ch	50
LY-3700	Maxi Tip	275mm	22ch	50
LY-3700/F	Maxi Tip (With blue Female connector)	315mm	22ch	50
LY-3700/RF/FF	Maxi Tip (with straight connector)	315mm	22ch	50
LY-3701	Maxi Tip (VC)	275mm	22ch	50
LY-3702	Maxi Tip (No side eyes)	275mm	22ch	50
LY-4700	Maxi Super Tip	275mm	22ch	50
LY-4700/PR	Maxi Super Tip (PR)	275mm	22ch	50
LY-4701	Maxi Super Tip (VC)	275mm	22ch	50
LY-3800	Maxi Plus Tip	280mm	28ch	50
LY-3801	Maxi Plus Tip (VC)	280mm	28ch	50
LY-3803	Maxi Plus Tip (VC) (No side eyes)	280mm	28ch	50



© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.