


# Link Yankauer Sets

## Product Data Sheet

<b>Product Codes</b>	See table over page
<b>Description</b>	Single use, Link Suction Yankauers, pre-connected to Pennine suction connecting tubing in 2m, 2.2m and 3metre lengths.
<b>Intended Use</b>	Surgical suction probe range of devices are hand operated, single use, and are used as a conduit to remove bodily fluids and secretions, surgical tissue debris, irrigation fluids and gasses. They are used in surgical sites and body cavities to facilitate observation and/or to clear an obstruction. To be used on all patient populations by clinicians and trained homecare users, in clinical setting, homecare and emergency settings.
<b>Duration of Use</b>	Normally intended for continuous use of less than 60 minutes, but will be patient specific depending on the duration of the surgical procedure.
<b>Product Information</b>	Available with or without vacuum control (VC) vent Short mini, Mini fine, Midi and Maxi Yankauers available with smooth rounded end or Super Tip. Tubing available in 2m, 2.2m and 3m lengths
<b>Materials</b>	Tube - PVC (DEHP Free) Funnel Connector - PVC (DEHP Free) Yankauer - PVC (DEHP Free) Latex Free
<b>Quality System</b>	ISO 13485
<b>EU Device Classification</b>	Class IIa
<b>Box Qty</b>	See table over page
<b>Packaging</b>	Sterile, double wrapped individual peel pouches
<b>Sterilised by</b>	Ethylene Oxide;
<b>Shelf Life</b>	5 years
<b>Storage &amp; Safety</b>	Toxicity: Manufactured from biocompatible materials. Storage: Store in a clean, cool, dry place away from direct sunlight



© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

# Link Yankauer Sets

## Product Data Sheet

### Sets with 2m & 2.2m tubing

Product Code	Tubing Diam/Length	Qty Per Box	Yankauer Type
LYS-5520	6mm / 2m	10	LY-3500 Mini Fine Tip
LYS-5620	6mm / 2m	10	LY-3600 Midi
LYS-5621	6mm / 2m	10	LY-3601 Midi VC
LYS-5720	6mm / 2m	10	LY-3700 Maxi
LYS-5721	6mm / 2m	10	LY-3701 Maxi VC
LYS-6620	7mm / 2.2m	10	LY-3600 Midi

### Sets with 3m tubing

Product Code	Tubing Diam/Length	Qty Per Box	Yankauer Type
LYS-5430	6mm / 3m	10	LY-4400 Short Mini
LYS-5530	6mm / 3m	10	LY-3500 Mini Fine Tip
LYS-5531	6mm / 3m	10	LY-3501 Mini Fine Tip VC
LYS-5630	6mm / 3m	10	LY-3600 Midi
LYS-5631	6mm / 3m	10	LY-3601 Midi VC
LYS-5730	6mm / 3m	10	LY-3700 Maxi
LYS-5731	6mm / 3m	10	LY-3701 Maxi VC
LYS-6830	7mm / 3m	5	LY-3800 Maxi Plus
LYS-7730	8mm/3m	5	LY-3700 Maxi



© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.