

Mucus Extractors

Product Data Sheet

Product Codes	See over page
Description	Mucus Extractor designed to obtain mucus specimen for micro-biological examination.
Intended Use	A single-use, covered plastic receptacle and associated tubing, containing no additives intended to be used for the aspiration, and/or collection of mucus/sputum specimens required for further examination. To be used on all patient populations, by clinician or suitably trained person, in a clinical setting.
Duration of Use	Normally intended for continuous use of less than 60 minutes.
Product Information	Fitted with 33cm suction catheter Collection Bottle Graduated to 20ml Supplied with spare screw cap Suction tube with standard funnel or vacuum control connector for attachment to suction source
Materials	Collection Bottle - Phillips K-resin Collection Bottle Cap - Polypropylene Tubing and Connector - PVC (DEHP Free) Latex Free
Quality System	ISO 13485
EU Device Classification	Class Is
Box Qty	See over page
Packaging	Sterile, individual peel pouch
Sterilised by	Ethylene Oxide;
Shelf Life	5 years
Shore Hardness of Connecting Tube	77a
Product CH	[See Overleaf]
Storage & Safety	Storage: Store in a clean, cool, dry place away from direct sunlight



© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.

Mucus Extractors

Product Data Sheet



Without Vacuum Control Connector

Product Code	Suction Tube Length	CH/Diam	Qty per box
ME-6004	16cm	14ch ø 4.66mm	100
ME-8108	25cm	8ch ø 2.66mm	75
ME-8110	25cm	10ch ø 3.33mm	75
ME-8112	25cm	12ch ø 4.00mm	75
ME-8114	25cm	14ch ø 4.66mm	75



With Vacuum Control Connector

Product Code	Suction Tube Length	CH/Diam	Qty per box
ME-9108	25cm	8ch ø 2.66mm	75
ME-9110	25cm	10ch ø 3.33mm	75
ME-9112	25cm	12ch ø 4.00mm	75
ME-9114	25cm	14ch ø 4.66mm	75
ME-9214 (extra long cannula 50cm)	25cm	14ch ø 4.66mm	75

© 2021 Pennine Healthcare or one of its subsidiaries. All rights reserved.

The information given in this sheet is not intended to be exhaustive.

The information contained in this sheet is liable to modification from time to time in the light of experience and our policy of continuous product development.