

QUFORA® LEG BAGS

Qufora® Leg bags are comfortable to wear, having a nice fit to the leg. The bag is almost invisible underneath the clothes when filled and the soft backing ensures high comfort, even in warm conditions.

Qufora® Leg bags are subject to strict quality requirements to meet the expectations of a user-friendly and durable bag.

PRODUCT FEATURES

- ▶ Five different sizes – 500 ml and 750 ml with a multi-chamber system as well as 350 ml, 600 ml, 750 ml and 800 ml without multi-chamber system.
- ▶ Sterile and non-sterile version. Most of the bags have an adjustable tube.
- ▶ Sterile Leg bags are available with needle-free outlet for urine sampling in four sizes.
- ▶ The universal connector makes it easy to connect the Leg bag to a catheter or a male external catheter.
- ▶ The discreet, soft and flexible tube ensures a high degree of freedom of movement.
- ▶ The non-return valve prevents the urine from running back in the tubing, thus preventing infection risk.
- ▶ The multi-chamber system ensures that the Leg bag is almost invisible under clothing, even when full.
- ▶ The back side has a skin-friendly 'non-woven' soft backing for increased comfort.
- ▶ Available with either the discrete Q-Flow valve, T-valve or Lever valve - all can be operated single-handed.
- ▶ The drainage valves are easy to operate and can be emptied in a hygienic and safe manner.
- ▶ The Leg bag can be connected directly to a Bed bag. This reduces the risk of infection in the event of rupture between catheter and Leg bag.
- ▶ The Leg bags are individually packaged.


Flexible tube with non-return valve



Holes for attaching velcro strap



Scale and comfortable soft backing on multi chamber bags



Q-flow or T-tap valve, operated with one hand



Connectors - available in two versions - with or without needle-free sampling





NON STERILE

MULTI CHAMBER - DRAINABLE - ADJUSTABLE TUBE / LOOSE CONNECTOR / CAP - SINGLE PACKED

REF	Volume	Tube	Valve	Soft backing	Qty/box
11121501	500 ml	50 cm	T-valve	Yes	10
11121502 *	500 ml	50 cm	Q-Flow	Yes	10
11121802 **	500 ml	10 cm	Q-Flow	Yes	10
11161502	500 ml	50 cm	Q-Flow	Yes	10
11321501	750 ml	50 cm	T-valve	Yes	10
11361502	750 ml	50 cm	Q-Flow	Yes	10

*with 1 pair of straps; ** no adjustable tube.

STANDARD - DRAINABLE - SINGLE PACKED

REF	Volume	Tube	Valve	Soft backing	Qty/box
11061901 *	350 ml	30 cm	Q-Flow	Yes	10
11221502 *	600 ml	50 cm	Q-Flow	Yes	25
11221902	600 ml	30 cm	Q-Flow	Yes	25
11223802 **	600 ml	10 cm	Q-Flow	No	25
11251502 *	600 ml	50 cm	T-valve	Yes	25
11251902	600 ml	30 cm	T-valve	Yes	25
11259802 **	600 ml	10 cm	T-valve	No	25
11451502	750 ml	50 cm	T-valve	No	25
11451803	750 ml	10 cm	T-valve	No	25
11421502 *	800 ml	50 cm	Lever	Yes	25
11451502 *	800 ml	50 cm	T-valve	Yes	25

*adjustable tube / loose connector / cap; **slanting tube.

COMBI - DRAINABLE - LOOSE CONNECTOR / CAP - SINGLE PACKED

REF	Volume	Tube	Valve	Soft backing	Qty/box
11751101 *	1.500 ml	90 cm	T-valve	Yes	10

*adjustable tube / loose connector / cap

STERILE

DUO – DRAINABLE - NEEDLE-FREE SAMPLE PORT - SINGLE PACKED

REF	Volume	Tube	Valve	Soft backing	Qty/box
12561001 *	1.000 ml	90 cm	Q-Flow	Yes	20

* extra needle-free sample port; "doghouse"; fixation string; hanger; odor filter; wide inlet tube

MULTI CHAMBER – DRAINABLE - NEEDLE-FREE SAMPLE PORT - SINGLE PACKED

REF	Volume	Tube	Valve	Soft backing	Qty/box
12061901 *	350 ml	30 cm	Q-Flow	Yes	10
12161502 *	500 ml	50 cm	Q-Flow	Yes	10
12161802	500 ml	10 cm	Q-Flow	Yes	10
12221502 *	600 ml	50 cm	Lever	Yes	10
12251902 *	600 ml	30 cm	T-valve	Yes	10
12361502 *	750 ml	50 cm	Q-Flow	Yes	10
12361802	750 ml	10 cm	Q-Flow	Yes	10

* adjustable tube / loose connector/cap; 12221502: no needle-free sample port; 12251902 and 12161802: wide inlet tube