

QUFORA® BED BAGS

Qufora® Bed bags are designed for hospitals, nursing homes and home care with focus on functionality and security.

Qufora® Bed bags are available in several versions and sizes to suit individual needs.

PRODUCT FEATURES

- ▶ A wide product assortment. Available in sterile or non-sterile, long or short tubes, drainable and non-drainable and with or without needle-free sampling.
- ▶ Has a universal connector which easily connects to the catheter or Leg bag.
- ▶ The flexible connection tube provides for safe drainage and allows the bag to be placed below bladder level.
- ▶ The non-return valve prevents the urine from running back in the tubing, thus preventing infection risk.
- ▶ Capacity: 750ml, 1.500 ml or 2.000 ml.
- ▶ Two different drainage valves - free-flow T-valves or large Lever valve. The drainage valves are easy to operate.
- ▶ The Bed bag can be emptied in a hygienic and safe way.
- ▶ Clear scale on the outside makes it easy to read the urine volume.
- ▶ The free-flow or large outlet ensures that clots freely can pass.
- ▶ Tube lengths from 10 cm - 130 cm.

“Eyes” for attaching to a bed bag stand



Connectors - available in two versions - with or without needle-free sampling



Both drainage valves are easy to operate with one hand




Non-drainable bag



“Doghouse”





NON STERILE

NON-DRAINABLE – TEAR MARK – SINGLE PACKED

REF	Volume	Tube	Valve	Single packed	Qty/box
13303802	750 ml	10 cm	-	Yes	25
13703412	1.500 ml	110 cm	-	Yes	25
13703802	1.500 ml	10 cm	-	Yes	25
13703902	1.500 ml	35 cm	-	Yes	25
13800101	2.000 ml	90 cm	-	No	10
13803112	2.000 ml	90 cm	-	Yes	25

DRAINABLE – SINGLE PACKED

REF	Volume	Tube	Valve	Single packed	Qty/box
13753101	1.500 ml	90 cm	T-valve	Yes	25
13850301	2.000 ml	130 cm	T-valve	Yes	10
13851202 *	2.000 ml	130 cm	Lever	Yes	10
13853101	2.000 ml	90 cm	T-valve	Yes	25
13853401	2.000 ml	110 cm	T-valve	Yes	25
13853402 **	2.000 ml	110 cm	T-valve	Yes	100
13853403 **	2.000 ml	110 cm	T-valve	Yes	100

* needle-free sample port; extra needle-free sample port; odor filter; hanger; sliding clamp; bed sheet clip; wide inlet tube

** wide inlet tube

STERILE

DRAINABLE – NEEDLE-FREE SAMPLE PORT – ODOR FILTER - SINGLE PACKED

REF	Volume	Tube	Valve	Single packed	Qty/box
14823202 – U4	2.000 ml	120 cm	Lever	Yes	25
14851201 – S4	2.000 ml	120 cm	Lever	Yes	25
14853202 – Q4	2.000 ml	120 cm	T-valve	Yes	25
14853203 *	2.000 ml	120 cm	T-valve	Yes	50
16823201 **	2.000 ml	130 cm	Lever	Yes	20

- U4/Q4: fold-away holder; bed sheet clip; tube clamp; wide inlet tube

- S4: "dog-house"; bed sheet clip; tube clamp; wide inlet tube

* no odor filter, bed sheet clip or tube clamp

** closed system: drip chamber; "doghouse"; hook; odor filter; fixation string; bed sheet clip