

## evercare<sup>®</sup> Medikit Forceps Helstead Mosquito

REF: 350152-1

Single-use sterile instruments are medical tools designed for one-time use during surgical procedures. Unlike traditional reusable instruments, which require cleaning, sterilization, and maintenance, single-use instruments are discarded after a single procedure. Single-use metal sterile instruments, such as forceps or scissors, serve a vital purpose in medical settings, primarily aiding in surgical procedures and various medical interventions. These instruments play a crucial role in modern healthcare settings, ensuring optimal sterility and reducing the risk of infections and cross-contamination. Single-use instruments are used on an individual patient during a single procedure and then disposed of, ensuring optimal hygiene and reducing the risk of cross-contamination between patients. Additionally, single-use metal sterile instruments offer convenience and cost-effectiveness for medical facilities, as they eliminate the need for expensive sterilization equipment and the associated maintenance costs. Healthcare providers can focus on patient care without the logistical challenges of instrument reprocessing, making these instruments an efficient solution for surgical procedures and medical interventions.

### Key Features

- ▶ sterile with each use,
- ▶ eliminating the need for extensive cleaning and sterilization processes,
- ▶ saving time,
- ▶ mitigating errors in the sterilization process,
- ▶ enhancing patient's safety.



REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
350152-1	evercare <sup>®</sup> Medikit Forceps Helstead Mosquito +non woven swab 5x5, White/Silver	13 cm	100 / 400

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-7
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR Class IIa
<b>Compliance with product standards:</b>	EN ISO 11135
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 20417, EN ISO 11607-1, EN ISO 15223-1

## Material

Product Part	Raw Material & Characteristics
100% stainless steel	stainless steel
non woven	viscose and polyester non woven 30 g/sqm

## Warnings & Recommendations

Warnings	
<b>Storage Recommendation</b>	Recommended storage in clean, dry space between 5 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, PP/PE film	200 x 120 x 10
<b>Dispenser Pack</b>	Cardboard	188 x 140 x 320
<b>Transport Carton</b>	TCF or ECF quality	385 x 286 x 335
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
<b>350152-1</b>	03666002021671	03666002021664	03666002021657

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.