

## Cotton

**100 g pleated cotton**  
**Width 12 cm**



- Highly absorbent 100% pure cotton for medical use
- Produced according to European Pharmacopoeia
- Refined and packed in Europe
- Appropriate for skin cleansing
- Not for wound cleaning or dressing directly on damaged skin
- Nonsterile
- Article number: BV0003-00

### **Packing**

Packed in resealable plastic bag  
100 g per package  
100 packages per transport box

### **Typical values**

pH value	4,4-8,2
Foreign fiber	occasional
Absorption cap.	>23g/g cotton
Sinking time	<10 sek
Micronaire	>3,5
Moisture	6-8 %
Ash content	<0,4%