

Cotton

500 g pleated cotton



- Highly absorbent 100% pure cotton for medical use
- Produced according to European Pharmacopoeia
- Refined and packed in Europe
- Appropriate for skin cleansing
- Not for wound cleaning or dressing directly on damaged skin
- Nonsterile

Packing

Packed in resealable plastic bag
500 g per package
20 packages per transport box

Typical values

pH value	4,4-8,2
Foreign fiber	occasional
Absorption cap.	>23g/g cotton
Sinking time	<10 sek
Micronaire	>3,5
Moisture	6-8 %
Ash content	<0,4%