


Cotton knitted gloves, Basic, natural white REF 222895S, 222895M, 222895L / S/M/L

Version: rev.1
Revision date: 06.06.2013

Last modified: 16.8.2013

| Sizes | REF | |
|--|--|--|
| S M L | 222895S 222895M 222895L |  |
| Colour: | natural white | |
| Overall length: | 30 cm | |
| <p>[NATURAL LATEX FREE]</p> | | |
| Raw material: | 100 % cotton | |
| REACH regulations: | The products covered by this data sheet do not contain any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA). | |
| Features: | Comfortable gloves made of 100 % cotton. Used e.g. in hand therapy (to assist absorption of creams and ointments, to prevent staining, to conceal skin disorders), also used as liners underneath waterproof gloves (e.g. protective gloves and household gloves). | |
| Pack count: | 40 pcs/bag, 50 bags (=2000 pcs)/transport carton | |
| Sales unit: | Transport carton | |
| Dispenser bag: | Polyethylene (PE) 330 x 140 mm | |
| Transport carton: | Corrugated cardboard S: 440 x 300 x 420 mm M: 460 x 300 x 420 mm L: 480 x 300 x 420 mm | |
| Compliance with legal requirements: | Medical Device Directive 93/42/EEC, class I | |
| Quality standards: | The product is manufactured in compliance with EN ISO 13485 | |
| Label information: | In accordance with EN980, EN ISO 15223-1 and EN1041 | |
| Storage recommendation: | Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone. | |
| Disposal recommendation: | Follow local instructions | |
| Risk if reused: | The product is single use. If reused, the performance of product may deteriorate. | |

Cotton knitted gloves, Basic, natural white REF 222895S, 222895M, 222895L / S/M/L

Page 2 (2)

Version: rev.1
Revision date: 06.06.2013

Last modified: 16.8.2013

| Barcodes: REF | Dispenser | Transport carton |
|----------------|---|--|
| 222895S | 6438129350838 222895S dispenser  6 438129 350838 | 6438129550832 222895S carton  6 438129 550832 |
| 222895M | 6438129350845 222895M dispenser  6 438129 350845 | 6438129550849 222895M carton  6 438129 550849 |
| 222895L | 6438129350852 222895L dispenser  6 438129 350852 | 6438129550856 222895L carton  6 438129 550856 |