



General description


- Made of thermosensitive PCV
- Safeguard to avoid tuck into the epipharynx
- Siliconised
- Latex free, phthalate free
- Single use

Design features

SIZE/INNER Ø ID [MM]	OUTER Ø OD [MM]	LENGTH L [MM]	REFERENCE NUMBER
2,5	3,5	81	41-2500
3,0	4,2	81	41-3000
3,5	4,8	91	41-3500
4,0	5,5	101	41-4000
4,5	6,0	111	41-4500
5,0	6,9	121	41-5000
5,5	7,5	131	41-5500
6,0	8,2	141	41-6000
6,5	8,7	141	41-6500
7,0	9,4	161	41-7000
7,5	10,2	161	41-7500
8,0	11,0	171	41-8000
8,5	11,7	171	41-8500
9,0	12,3	181	41-9000



PART NO.	PART NAME	MATERIAL
1	Tube body	PVC
2	Flange	PVC



Delivery

Standard delivery :

- Package dimension:
399 x 399 x 360 mm
- Quantity – 100 pcs.
- Approximate weight of the package:
2,50 – 3,50 kg

Custom delivery

Minimum order quantity – 10 pcs.