


# evercare<sup>®</sup> Hand/foot drape 230 x 300 cm, aperture 3,5 cm

REF: 1324-02

evercare<sup>®</sup> Hand/foot drape, sterile with extra absorbent patch, tube holders and elastic aperture for procedures in knees and feet.

## Key Features

- ▶ Absorbent patch
- ▶ Tube holders
- ▶ Elastic aperture
- ▶ Complies to EN 13795-1 High performance critical product area requirements



230 x 300 cm

REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1324-02	Hand/foot drape, sterile, Blue	230 x 300 cm, elastic Ø 3.5 cm	8 / 16

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



**REACH regulation:** The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-01  
EN ISO 10993-5  
EN ISO 10993-7  
EN ISO 10993-10

**Compliance with legal requirements:** Medical Device Regulation (EU) 2017/745, MDR class Is

**Compliance with product standards:** EN13795-1

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
Laminate	PP nonwoven, blue / PE film, blue / PP nonwoven, white
Absorbent patch	PP nonwoven, dark blue
Elastic plastic film	SEBS-based thermoplastic elastomer
Tube holders	PP blue

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
Single Pack	Medical grade paper, PE/PET film	420 x 603
Dispenser Pack	TCF or ECF quality	485 x 285 x 287
Transport Carton	TCF or ECF quality	500 x 310 x 590
Packing Tape	Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1324-02	06438129119084	06438129319088	06438129519082

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.