


evercare[®] Stockinette, large, 30 x 120 cm

REF: 1395-01

evercare[®] Stockinette, sterile is an easy to drape two-layer stockinette to cover limbs before draping, particularly during orthopedic surgery. Outer layer elastic membrane, inner layer soft and comfortable.



30 x 120 cm

Key Features


- ▶ Easy to drape
- ▶ Inner layer soft and comfortable

REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1395-01	Stockinette, sterile, Blue	Large, 30 x 120 cm	10 / 40

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-01
EN ISO 10993-5
EN ISO 10993-7
EN ISO 10993-10

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR class Is

Compliance with product standards: EN ISO 11135

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Top	Synthetic rubber, blue
Lining	Elastic gauze, PES

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
Storage Recommendation	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PE/PET film	250 x 390
Dispenser Pack	TCF or ECF quality	285 x 242 x 287
Transport Carton	TCF or ECF quality	500 x 310 x 590
Packing Tape	Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
1395-01	06438129118896	06438129318890	06438129518894

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.