


### Elastic stockinettes

#### REF 1395-01, 1396-02

REF	description	size																									
1395-01	Stockinette	30 x 120 cm, L																									
1396-01	Stockinette	22 x 75 cm, S																									
CE0434 The product complies with legal requirements for medical devices EO sterilized according to valid version of ISO 11135 standard <b>[ NATURAL LATEX FREE ]</b>																											
<b>Raw materials</b> Top Lining			Synthetic rubber, blue Elastic gauze, PES																								
<b>REACH regulations:</b>			The products covered by this data sheet do not contain carcinogenic, mutagenic or toxic substances, including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).																								
<b>Features:</b>			Easy to drape two-layer stockinet's to isolate limbs particularly during orthopedic surgery. Outer layer impervious, inner layer soft and comfortable.																								
<b>Shelf life:</b>			5 years																								
<b>Pack count, dispenser/carton measurements</b>			<table border="1"> <thead> <tr> <th>REF</th> <th>pcs/ unit pack</th> <th>pcs/ dis- penser</th> <th>pcs/ transp. ctn</th> <th>meas. unit pack mm</th> <th>meas. disp. mm</th> <th>meas. transp. ctn mm</th> </tr> </thead> <tbody> <tr> <td>1395-01</td> <td>1</td> <td>10</td> <td>40</td> <td>250x390</td> <td>285x242x287</td> <td>500x300x590</td> </tr> <tr> <td>1396-01</td> <td>1</td> <td>20</td> <td>80</td> <td>150x289</td> <td>265x143x246</td> <td>545x300x246</td> </tr> </tbody> </table>				REF	pcs/ unit pack	pcs/ dis- penser	pcs/ transp. ctn	meas. unit pack mm	meas. disp. mm	meas. transp. ctn mm	1395-01	1	10	40	250x390	285x242x287	500x300x590	1396-01	1	20	80	150x289	265x143x246	545x300x246
REF	pcs/ unit pack	pcs/ dis- penser	pcs/ transp. ctn	meas. unit pack mm	meas. disp. mm	meas. transp. ctn mm																					
1395-01	1	10	40	250x390	285x242x287	500x300x590																					
1396-01	1	20	80	150x289	265x143x246	545x300x246																					
<b>Sales unit:</b>			Transport carton																								
<b>Unit pack:</b>			Medical grade paper/PE film pouch																								
<b>Dispenser box</b>			Cardboard, totally chlorine free																								
<b>Transport carton:</b>			Cardboard, totally chlorine free																								
<b>Packing tape:</b>			Transparent, PVC-free																								
<b>Compliance with legal requirements:</b>			Medical Devices Directive 93/42/EEC, Class Is																								
<b>Compliance with product standards:</b>			EN13795																								
<b>Quality standards:</b>			The product is manufactured in compliance with, ISO13485																								
<b>Label information:</b>			in accordance with EN ISO 15223-1 and EN1041																								
<b>Storage recommendation:</b>			Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone																								
<b>Disposal recommendation</b>			Follow local instructions																								
<b>Risk if reused:</b>			The product is single use. If reused performance of product may deteriorate, cross contamination may occur																								
<b>Barcodes:</b> Unit pack		Dispenser	Transport carton																								

### Elastic stockinettes REF 1395-01, 1396-02

Version: 20.10.2015, rev 1

Last modified: 26.10.2015

6438129118896 1395-01 single  6 438129 118896	6438129318890 1395-01 dispenser  6 438129 318890	6438129518894 1395-01 carton  6 438129 518894
6438129119466 1396-01 single  6 438129 119466	6438129319460 1396-01 dispenser  6 438129 319460	6438129519464 1396-01 carton  6 438129 519464