


# evercare<sup>®</sup> Adhesive aperture drape, variable, 50 x 75 cm

REF: 1512-01

evercare<sup>®</sup> Adhesive aperture drapes, variable, sterile with integrated edge fixation provides easy and quick draping for a variety of procedures.

## Key Features

- ▶ Easy and quick draping
- ▶ Complies to EN 13795-1 High performance critical product area requirements




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1512-01	Adhesive aperture drapes, variable, sterile, Blue	50 x 75 cm, adhesive edge 5 cm	35 / 140

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



**REACH regulation:** The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-01  
EN ISO 10993-5  
EN ISO 10993-7  
EN ISO 10993-10

**Compliance with legal requirements:** Medical Device Regulation (EU) 2017/745, MDR class Is

**Compliance with product standards:** EN13795-1

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
Laminate	PP nonwoven , blue / PE-film, blue

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, PE/PET film	250 x 390
<b>Dispenser Pack</b>	TCF or ECF quality	285 x 242 x 287
<b>Transport Carton</b>	TCF or ECF quality	500 x 310 x 590
<b>Packing Tape</b>	Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
<b>1512-01</b>	06438129116229	06438129316223	06438129516227

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.