


# evercare<sup>®</sup> Adhesive aperture drape 75 x 90 cm, aperture 6 x 8 cm

REF: 1515-01

evercare<sup>®</sup> Adhesive aperture drape, sterile with adhesive aperture provides easy and quick draping for a variety of procedures.

## Key Features

- ▶ Easy and quick draping
- ▶ Complies to EN 13795-1 High performance critical product area requirements




75 x 90 cm  
 ø 6 x 8 cm

REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
1515-01	Adhesive aperture drape, sterile, Blue	75 x 90 cm, ø 6 x 8 cm	25 / 100

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



**REACH regulation:** The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-01  
EN ISO 10993-5  
EN ISO 10993-7  
EN ISO 10993-10

**Compliance with legal requirements:** Medical Device Regulation (EU) 2017/745, MDR class Is

**Compliance with product standards:** EN13795-1

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
Laminate	PP nonwoven, blue / PE-film, blue

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur. Do not use if package is damaged.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, Transparent	150 x 290
<b>Dispenser Pack</b>	TCF or ECF quality	265 x 143 x 234
<b>Transport Carton</b>	TCF or ECF quality	545 x 310 x 246
<b>Packing Tape</b>	Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
<b>1515-01</b>	06438129110647	06438129310641	06438129510645

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.