


1870 Fluid collection pouches

Page 1 (2)

Version: 27.02.2014, rev. 1

Last modified: 27.2.2014


REF	Description	Size
		cm
1870	Fluid collection pouch, with adhesive edge	30x40

<p>CE0434 The product complies with legal requirements for medical devices EO sterilized according to valid version of ISO 11135 standard [NATURAL LATEX FREE]</p>	 <p>30 x 40 cm</p>
<p>Contents</p> <p>Main raw material Sieve Outlet valve & capp Support wire</p>	<p>Raw material</p> <p>100 % PEVA, transparent PE/PET, white 100 % PE PE-coated metal wire</p>
<p>REACH regulations:</p>	<p>The products covered by this data sheet specification do not contain carcinogenic, mutagenic or toxic substances, including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).</p>
<p>Features:</p>	<p>Pouches for collection of substances such as blood, bodily fluids, rinsing fluids, body tissue from operation site etc</p>
<p>Shelf life:</p>	<p>5 years</p>
<p>Pack count, carton weight:</p>	<p>1870 1 pc/bag, 55 pcs/dispenser, 220 (4 x 55)/transport carton</p>
<p>Sales unit:</p>	<p>Transport carton</p>
<p>Unit pack:</p>	<p>Medical grade paper/PE film pouch</p>
<p>Dispenser box</p>	<p>1870 Cardboard, totally chlorine free, size 143x234x265 mm</p>
<p>Transport carton:</p>	<p>1870 Cardboard, totally chlorine free, size 300x246x545 mm</p>
<p>Packing tape:</p>	<p>Transparent, PVC-free</p>
<p>Compliance with legal requirements:</p>	<p>Medical Devices Directive 93/42/EEC, Class Is</p>
<p>Compliance with product standards:</p>	<p>EN13795</p>
<p>Quality standards:</p>	<p>The product is manufactured in compliance with, EN ISO 13485 in accordance with EN ISO 15223-1 EN 980 and EN1041</p>
<p>Label information:</p>	
<p>Storage recommendation:</p>	<p>Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone</p>
<p>Disposal recommendation</p>	<p>Follow local instructions</p>
<p>Risk if reused:</p>	<p>The product is single use. If reused performance of product may deteriorate, cross contamination may occur</p>

1870 Fluid collection pouches

Version: 27.02.2014, rev. 1

Last modified: 27.2.2014

Barcodes: Unit pack	Dispenser	Transport carton
6438129111736 1870 single  6 438129 111736	6438129311730 1870 dispenser  6 438129 311730	6438129511734 1870 carton  6 438129 511734