

# Product Data Sheet

Staff clothing

Surgical gowns evercare<sup>®</sup> XP,  
Standard

REF 2504 - 2510



Page 1 (2)

Version: rev.1,15.05.2012

Printed: 5.7.2012

REF 2504	120 cm, M-short, without hand towels	
REF 2505	130 cm, M, without hand towels	
REF 2507	140 cm, L, without hand towels	
REF 2509	150 cm, XL, without hand towels	
REF 2510	160 cm, XXL, without hand towels	
CE0434	The product complies with legal requirements for medical devices	
EO sterilized according to valid version of ISO 11135-1:2007 standard		
<div style="border: 1px solid blue; padding: 2px; display: inline-block;"><b>NATURAL LATEX FREE</b></div>		
<b>Product part</b> Main raw material Cuffs Neck ribbon Velcro band Belt card Sewing thread Crepe wrapping Unit pack Dispenser box Transport carton	<b>Raw material</b> Blue SMS nonwoven, multilayer White rib knitted fabric White spunlace nonwoven Woven ribbon with hook or loop White White. Natural white paper Paper/plastic White cardboard Brown cardboard	<b>Composition</b> PP PES PES PA Pulp (ECF) PES Pulp (ECF) Medical grade paper (TCF)/LLDPE Pulp (TCF) Pulp (TCF)
<b>REACH regulations:</b>	The product covered by this data sheet does not contain any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).	
<b>Features:</b>	Surgical gown with overlapping back, the "universal" gown for dry operations	
<b>Shelf life:</b>	5 years	
<b>Pack count,</b>	REF 2504	1/21/42
	REF 2505	1/20/40
	REF 2507	1/20/40
	REF 2509	1/18/36
	REF 2510	1/18/36
<b>Sales unit</b>	Transport carton	
<b>Unit pack</b>	Medical grade paper (TCF)/LLDPE pouch 250 x 388 mm	
<b>Dispenser box</b>	TCF cellulose, 381 x 279 x 320 mm	
<b>Transport carton:</b>	TCF cellulose, 574 x 391 x 330 mm	
<b>Packing tape:</b>	PVC-free	
<b>Compliance with legal requirements:</b>	Medical Devices Directive 93/42/EEC, Class I	
<b>Compliance with product standards:</b>	EN13795	
<b>Quality standards:</b>	The product is manufactured in compliance with ISO13485	
<b>Label information:</b>	in accordance with EN980 and EN1041	
<b>Storage recommendation:</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight , other intensive light sources and ozone	
<b>Disposal recommendation</b>	Energy recovery or mixed waste	

# Product Data Sheet

Staff clothing

Surgical gowns evercare<sup>®</sup> XP,  
Standard

REF 2504 - 2510



Page 2 (2)

Version: rev.1,15.05.2012

Printed: 5.7.2012

REF	Unit pack	Dispenser	Transport carton
<b>REF 2504</b>	6438129112139 2504 single  6 438129 112139	6438129312133 2504 dispenser  6 438129 312133	6438129512137 2504 carton  6 438129 512137
<b>REF 2505</b>	6438129112146 2505 single  6 438129 112146	6438129312140 2505 dispenser  6 438129 312140	6438129512144 2505 carton  6 438129 512144
<b>REF 2507</b>	6438129112153 2507 single  6 438129 112153	6438129312157 2507 dispenser  6 438129 312157	6438129512151 2507 carton  6 438129 512151
<b>REF 2509</b>	6438129112160 2509 single  6 438129 112160	6438129312164 2509 dispenser  6 438129 312164	6438129512168 2509 carton  6 438129 512168
<b>REF 2510</b>	6438129112177 2510 single  6 438129 112177	6438129312171 2510 dispenser  6 438129 312171	6438129512175 2510 carton  6 438129 512175