

## evercare<sup>®</sup> Protective gown Special

REF: 2720

evercare<sup>®</sup> Protective gown Special is a comfortable water repellent protective gown for different areas of usage within healthcare. Made in a dark blue nonwoven material, with white comfortable cuffs and tie band in neck and waist to fit all sizes. The product is for single-use and non-sterile.

### Key Features

- ▶ Comfortable cuffs
- ▶ Tie band in neck and waist
- ▶ Water repellent
- ▶ 125 cm length









REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2720	Protective gown Special, 5-pack, Darkblue	L, length 125 cm	5 / 25

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices.

REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-05 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDD class I
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
Body, sleeve and belt	Spunbond PP nonwoven, dark blue
Cuff	PES, white
Neck ribbon	PES Spunlace nonwoven, white
Sewing thread	PES, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	PE bag	
<b>Transport Carton</b>	Carton, TCF or ECF quality	380 x 318 x 298
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
2720	N/A	N/A	6438129512601

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.