

# Surgical hood James

REF: 2773-01, 2774-01

**Surgical hood James** is a light version helmet without shoulder coverings. The surgical hood is made in a soft and breathable viscose nonwoven material (100% Fossile-free). The product is designed for staff who cannot use hoods covering the shoulders.

The use of headwear is proven to reduce the risk of particles falling into the surgical wound.

## Key Features:

- ▶ Surgical hood without shoulder protection
- ▶ Closure with strap underneath the chin, tied in the back
- ▶ Soft and breathable material, 100% fossile-free
- ▶ Comes in 2 colors: green and blue
- ▶ Universal size



REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
2773-01	Surgical hood James, green	Universal	100/600
2774-01	Surgical hood James, blue	Universal	100/600

**Sales Unit:** Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices



<b>REACH regulation:</b>	The product(s) covered by this data sheet does(do) not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Devices Directive 93/42/EEC, Class I
<b>Compliance with product standards:</b>	N/A
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristic
Body	Viscose nonwoven, green or blue
Sewing threads	Nylon, white

## Warnings and Recommendations



<b>Warnings:</b>	The product is single use. If reused performance of product may deteriorate, cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sun light and other intensive light sources.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 Years

## Packing Information

Package	Material	Size ( L x W x H ) mm
Single Pack	N/A	N/A

<b>Dispenser Pack</b>	Carton, TCF quality	145 x 185 x 300
<b>Transport Carton</b>	Carton, TCF quality	455 x 310 x 385
<b>Packing Tape</b>	Transparent, PVC-free	N/A

### Barcodes

REF	Single Pack	Dispenser	Transport Carton
2773-01	N/A	 6 438129 319651	 6 438129 519655
2774-01	N/A	 6 438129 319668	 6 438129 519662



Manufactured in China for  
**OneMed Group Oy**  
 Address: Metsäläntie 20  
 Postal Address: P.O. Box 10  
 FI-00321 Helsinki, Finland  
 Phone +358 20 786 6800  
 www.onemed.com  
 Domicile Helsinki  
 Business ID 2039640-1