

## evercare<sup>®</sup> Elastic support bandage

REF: 4251, 4252, 4253, 4254, 4255

evercare<sup>®</sup> Elastic support bandage is stretchable with soft and porous structure, have high permeability and is comfortable to wear. Thanks to woven-edge weaving technique and careful selection of raw material, the bandages may be used for support, fixation of most dressings, especially for joints and round body parts, and small medical devices to the patient's body. The bandage stays stable once applied and its layers do not slip. Package designed for easy product identification.

### Key Features

- ▶ Fixation of wound dressings, especially on joints and round body parts
- ▶ For support, e.g. after plaster removal and odema treatment
- ▶ Resistant to ointments
- ▶ Bandage with woven edges
- ▶ Soft and comfortable
- ▶ Light support bandage
- ▶ 60% stretchability
- ▶ 65% cotton and 35% polyamide




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4251	Elastic support bandage, cotton and polyamide	6 cm x 4 m	10 / 280
4252	Elastic support bandage, cotton and polyamide	8 cm x 4 m	10 / 240
4253	Elastic support bandage, cotton and polyamide	10 cm x 4 m	10 / 200
4254	Elastic support bandage, cotton and polyamide	12 cm x 4 m	10 / 160
4255	Elastic support bandage, cotton and polyamide	15 cm x 4 m	10 / 120

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-5 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Compliance with product standards:</b>	DIN 61632, DIN 61634
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN ISO 20417

## Material

Product Part	Raw Material & Characteristics
<b>Bandage</b>	65 % cotton and 35 % polyamide, natural white, Unstretched length 2.5 m, Stretched length 4 m, Elongation approx. 50%, Elasticity approx. 85%

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	5 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	4251: 205 x 85 x 65 4252: 205 x 85 x 85 4253: 205 x 85 x 105 4254: 205 x 85 x 125 4255: 205 x 85 x 155
<b>Transport Carton</b>	Carton, TCF or ECF quality	4251: 475 x 360 x 225 4252: 530 x 360 x 225 4253: 545 x 360 x 225 4254: 520 x 360 x 225 4255: 485 x 360 x 225
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4251	N/A	06438129314878	06438129514872
4252	N/A	06438129314885	06438129514889
4253	N/A	06438129314892	06438129514896
4254	N/A	06438129314908	06438129514902
4255	N/A	06438129314915	06438129514919

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.