

## evercare<sup>®</sup> Self-adhesive fabric

REF: 4555, 4560, 4565, 4570, 4575

evercare<sup>®</sup> Self-adhesive fabric is used to support fixation of dressings, swabs, catheters, tubes and cannulae. Skin-friendly adhesive, stretchable and conforms well to body contours. Allows movement when in place. Easy to remove release paper premeasured in 1 cm intervals. Breathable, pleasant to wear and moisture repellent.

### Key Features

- ▶ Ideal for areas such as neck, elbow and knee
- ▶ Perforated split for easy removal and application
- ▶ Barrier for blood and fluid strike through




REF	DESCRIPTION	SIZE	PACK COUNT dispenser/carton
4555	Self-adhesive fabric roll, White	2,5 cm x 10 m	1 / 120
4560	Self-adhesive fabric roll, White	5 cm x 10 m	1 / 60
4565	Self-adhesive fabric roll, White	10 cm x 10 m	1 / 36
4570	Self-adhesive fabric roll, White	15 cm x 10 m	1 / 24
4575	Self-adhesive fabric roll, White	20 cm x 10 m	1 / 18

Sales Unit: Transport Carton

# Technical Summary

## Regulations

The product complies with legal requirements for medical devices.



<b>REACH regulation:</b>	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
<b>Biological evaluation:</b>	EN ISO 10993-05 EN ISO 10993-10
<b>Compliance with legal requirements:</b>	Medical Device Regulation (EU) 2017/745, MDR class I
<b>Quality standards:</b>	EN ISO 13485
<b>Label information and packaging compliance:</b>	EN ISO 15223-1, EN 1041

## Material

Product Part	Raw Material & Characteristics
<b>Carrier</b>	100% polyester, nonwoven spunlance, water repellent
<b>Coating</b>	Pressure sensitive glue, thermosetting acrylic solution polymer, latex free
<b>Release paper</b>	Coated paper, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space in ambient temperature. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	3 Years

## Packing information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	N/A	N/A
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	4555: 180 x 87 x 30 4560: 30 x 83 x 83 4565: 55 x 83 x 83 4570: 105 x 83 x 83 4575: 155 x 83 x 83
<b>Transport Carton</b>	Carton, TCF or ECF quality	4555: 455 x 375 x 175 4560: 360 x 330 x 280 4565: 360 x 300 x 280 4570: 360 x 330 x 280 4575: 360 x 330 x 280
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes

REF	Single Pack	Dispenser	Transport Carton
<b>4555</b>	N/A	6438129313833	6438129513837
<b>4560</b>	N/A	6438129313840	6438129513844
<b>4565</b>	N/A	6438129313857	6438129513851
<b>4570</b>	N/A	6438129313864	6438129513868
<b>4575</b>	N/A	6438129313871	6438129513875

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.