

## evercare<sup>®</sup> Adhesive surgical dressing, sterile

REF: 4605, 4610, 4615, 4620, 4625, 4630, 4635

evercare<sup>®</sup> Adhesive surgical dressing, sterile is used for surgical wound treatment and exudate absorption. The absorbent pad is finished by semi-permeable net, allowing to absorb significant quantity of exudate coming from the wound while performing non-adherent properties at the same time. High retention rate allows to keep exudate inside absorbent pad even under external pressure. The dressing is air- and moisture vapor permeable. It can easily be applied on conical and round part of the patient's body.

### Key Features

- ▶ Soft and flexible
- ▶ High breathability
- ▶ For low to moderately exuding wounds
- ▶ Acrylic based adhesive
- ▶ Easy-to-use unit pouch
- ▶ Water repellent non-woven material




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4605	Adhesive surgical dressing with pad, sterile	5x7 cm; pad 2,5x4 cm	50 / 2400
4610	Adhesive surgical dressing with pad, sterile	10x10 cm; pad 5x5 cm	50 / 1200
4615	Adhesive surgical dressing with pad, sterile	10x15 cm; pad 5x11 cm	50 / 600
4620	Adhesive surgical dressing with pad, sterile	10x20 cm; pad 5x15,5 cm	30 / 360
4625	Adhesive surgical dressing with pad, sterile	10x25 cm; pad 5x20 cm	30 / 360
4630	Adhesive surgical dressing with pad, sterile	10x30 cm; pad 5x25 cm	30 / 360
4635	Adhesive surgical dressing with pad, sterile	10x35 cm; pad 5x30 cm	30 / 360

Sales Unit: Transport Carton

## Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



**REACH regulation:** The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

**Biological evaluation:** EN ISO 10993-5  
EN ISO 10993-7  
EN ISO 10993-10  
EN ISO 10993-11

**Compliance with legal requirements:** Medical Device Regulation (EU) 2017/745, MDR class IIa

**Compliance with product standards:** EN ISO 11135, EN ISO 11737-1, EN 13726-2

**Quality standards:** EN ISO 13485

**Label information and packaging compliance:** EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

## Material

Product Part	Raw Material & Characteristics
Carrier	100% polyester, nonwoven spunlance, water repellent, white
Coating	Pressure sensitive glue, thermosetting acrylic solution polymer, latex free
Absorbent pad	70 % viscose, 30 % polyester and non-adherent porous PE film, white
Release paper	Coated paper, white

## Warnings & Recommendations

<b>Warnings</b>	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
<b>Storage Recommendation</b>	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
<b>Disposal Recommendation</b>	Dispose in accordance with local regulations.
<b>Shelf Life</b>	3 years

## Packing Information

Package	Material	Size ( L x W x H ) mm
<b>Single Pack</b>	Medical grade paper, PP/PE film	4605: 72 x 113
		4610: 135 x 160
		4615: 135 x 220
		4620: 135 x 262
		4625: 135 x 313
		4630: 135 x 361
		4635: 135 x 413
<b>Dispenser Pack</b>	Carton, TCF or ECF quality	4605: 74 x 95 x 115
		4610: 95 x 140 x 165
		4615: 95 x 140 x 225
		4620: 65 x 140 x 265
		4625: 65 x 140 x 316
		4630: 65 x 140 x 365
		4635: 65 x 140 x 415
<b>Transport Carton</b>	Carton, TCF or ECF quality	4605: 320 x 600 x 260
		4610: 300 x 600 x 360
		4615: 300 x 600 x 240
		4620: 285 x 440 x 285
		4625: 285 x 440 x 335
		4630: 285 x 440 x 385
		4635: 285 x 440 x 435
<b>Packing Tape</b>	PVC-free, Transparent	N/A

## Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4605	06438129113891	06438129313895	06438129513899
4610	06438129113907	06438129313901	06438129513905
4615	06438129113914	06438129313918	06438129513912
4620	06438129113921	06438129313925	06438129513929
4625	06438129113938	06438129313932	06438129513936
4630	06438129113945	06438129313949	06438129513943
4635	06438129113952	06438129313956	06438129513950

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.