

evercare[®] Adhesive film dressing, sterile

REF: 4701, 4705, 4710, 4715, 4720, 4725

evercare[®] Adhesive film dressing, sterile is designed for wound treatment or cannula fixation. Adhesive film dressing maintains optimal moisture level of treated area while protecting against dust, water and bacteria at the same time. Polyurethane(PU) semi-permeable foil used for its composition allows the patient to take a shower reducing the risk of wound contamination. The dressing's transparency allows visual inspection of the wound without removal. The dressing design ensures painless removal.

Key Features

- ▶ Soft and flexible
- ▶ High breathability
- ▶ Acrylic based adhesive
- ▶ Easy-to-use unit pouch
- ▶ Effective barrier to exudate, water and bacteria
- ▶ Primary dressing for clean wounds
- ▶ Secondary dressing for fixation




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4701	Adhesive film dressing, sterile, Transparent	4cm x 5cm	100 / 2000
4705	Adhesive film dressing, sterile	6x7 cm	50 / 200
4710	Adhesive film dressing, sterile	10x12 cm	25 / 100
4715	Adhesive film dressing, sterile	10x25 cm	25 / 100
4720	Adhesive film dressing, sterile	15x20 cm	25 / 100
4725	Adhesive film dressing, sterile, Transparent	20x30 cm	10 / 100

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation: The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).

Biological evaluation: EN ISO 10993-5
EN ISO 10993-7
EN ISO 10993-10
EN ISO 10993-11

Compliance with legal requirements: Medical Device Regulation (EU) 2017/745, MDR class IIa

Compliance with product standards: EN ISO 11135, EN ISO 11737-1, EN 13726-2, EN 13726-4

Other tests performed: 0

Quality standards: EN ISO 13485

Label information and packaging compliance: EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Carrier	Transparent PU-film, breathable
Coating	Polyacrylic adhesive
Release paper	Medical grade coated paper

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	3 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	4701: 91 x 154 4705: 66 x 113 4710: 127 x 214 4715: 313 x 135 4720: 186 x 313 4725: 232 x 380
Dispenser Pack	Carton, TCF or ECF quality	4701: 70 x 50 x 115 4705: 95 x 30 x 160 4710: 130 x 20 x 220 4715: 140 x 20 x 315 4720: 190 x 20 x 320 4725: 240 x 15 x 382
Transport Carton	Carton, TCF or ECF quality	4701: 280 x 155 x 250 4705: 175 x 140 x 115 4710: 235 x 100 x 150 4715: 330 x 100 x 160 4720: 335 x 100 x 210 4725: 395 x 180 x 260
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4701	06438129124330	06438129324334	06438129524338
4705	06438129114164	06438129314168	06438129514162
4710	06438129114171	06438129314175	06438129514179
4715	06438129114188	06438129314182	06438129514186
4720	06438129114195	06438129314199	06438129514193
4725	06438129124347	06438129324341	06438129524345

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.