

evercare[®] Absorbent dressing, sterile

REF: 4781-01, 4782-01, 4783-01, 4784-01, 4785-01, 4786-01, 4787-01, 4788-01, 4789-01

evercare[®] Absorbent dressing, sterile is made of 4-ply: first a soft non-woven contact layer non-adherent to the wound, secondly an absorbent core made of cellulose fluff, thirdly a cellulose-tissue to encapsulate absorption pad, and finally an outer water-repellent non-woven layer which also protect the wound against external contamination. This design ensure high absorption capacity of the dressing while it retains exudate even under pressure.

Key Features

- ▶ High absorption and retention capacity of exudate
- ▶ Protects against external contamination
- ▶ Soft and flexible, can fit all body parts
- ▶ Non-adherent to the wound surface
- ▶ Steralized using EO
- ▶ Easy-to-use unit pouch




REF	DESCRIPTION	SIZE	Pack Count dispenser/carton
4781-01	Absorbent dressing, sterile	10x10 cm	50 / 600
4782-01	Absorbent dressing, sterile	10x20 cm	35 / 420
4783-01	Absorbent dressing, sterile	10x30 cm	25 / 300
4784-01	Absorbent dressing, sterile	15x20 cm	25 / 300
4785-01	Absorbent dressing, sterile	15x25 cm	25 / 300
4786-01	Absorbent dressing, sterile	20x20 cm	25 / 300
4787-01	Absorbent dressing, sterile	20x25 cm	20 / 240
4788-01	Absorbent dressing, sterile	20x30 cm	20 / 240
4789-01	Absorbent dressing, sterile	20x40 cm	10 / 120

Sales Unit: Transport Carton

Regulations

The product complies with legal requirements for medical devices. EO sterilized according to valid version of ISO 11135 standard.



REACH regulation:	The products covered by this data sheet do not contain substances classified as carcinogenic, mutagenic or toxic for reproduction (CMR), including any of the Substances of Very High Concern (SVHC) as listed in the latest available version of the Candidate List published by European Chemicals Agency (ECHA).
Biological evaluation:	EN ISO 10993-5 EN ISO 10993-7 EN ISO 10993-10
Compliance with legal requirements:	Medical Device Regulation (EU) 2017/745, MDR class Is
Compliance with product standards:	EN ISO 11135:2014/A1:2019
Quality standards:	EN ISO 13485
Label information and packaging compliance:	EN ISO 15223-1, EN ISO 20417, EN ISO 11607-1

Material

Product Part	Raw Material & Characteristics
Inner layer	PP nonwoven, white(this side against the wound)
Absorbent part	Wood pulp cotton encapsulated in fiber paper, white
Outer layer	PP nonwoven, blue, water repellent

Warnings & Recommendations

Warnings	The product is single use. If reused, the performance of product may deteriorate and cross contamination may occur.
Storage Recommendation	Recommended storage in clean, dry space between 15 - 25 °C. Products should be protected from direct sunlight, other intensive light sources and ozone.
Disposal Recommendation	Dispose in accordance with local regulations.
Shelf Life	5 years

Packing Information

Package	Material	Size (L x W x H) mm
Single Pack	Medical grade paper, PP/PE film	4781-01: 138 x 122
		4782-01: 249 x 132
		4783-01: 138 x 122
		4784-01: 194 x 133
		4785-01: 194 x 133
		4786-01: 249 x 132
		4787-01: 249 x 132
		4788-01: 249 x 132
		4789-01: 249 x 132
Dispenser Pack	Carton, TCF or ECF quality	4781-01: 140 x 120 x 275
		4782-01: 250 x 130 x 195
		4783-01: 140 x 120 x 340
		4784-01: 195 x 130 x 225
		4785-01: 195 x 130 x 340
		4786-01: 250 x 130 x 225
		4787-01: 250 x 130 x 270
		4788-01: 250 x 130 x 270
		4789-01: 250 x 130 x 270
Transport Carton	Carton, TCF or ECF quality	4781-01: 575 x 370 x 285
		4782-01: 765 x 270 x 400
		4783-01: 575 x 370 x 350
		4784-01: 600 x 270 x 460
		4785-01: 600 x 270 x 690
		4786-01: 765 x 270 x 460
		4787-01: 765 x 270 x 550
		4788-01: 765 x 270 x 550
		4789-01: 765 x 270 x 350
Packing Tape	PVC-free, Transparent	N/A

Barcodes / UDI-DI

REF	Single Pack	Dispenser	Transport Carton
4781-01	06438129121056	06438129321050	06438129521054
4782-01	06438129121063	06438129321067	06438129521061
4783-01	06438129121070	06438129321074	06438129521078
4784-01	06438129121087	06438129321081	06438129521085
4785-01	06438129121094	06438129321098	06438129521092
4786-01	06438129121100	06438129321104	06438129521108
4787-01	06438129121117	06438129321111	06438129521115
4788-01	06438129121124	06438129321128	06438129521122
4789-01	06438129121131	06438129321135	06438129521139

Evercare Medical AB is a partner that enables better lives for patients, facilitates everyday life for healthcare employees and contributes to lower care costs. We strive to reduce the environmental impact and, together with our suppliers, to improve working conditions throughout the value chain.